

# Analysis of Total Organic Nitrogen

File: TON SOP ver1.doc  
Last Modified: 12/13/2005 5:47:00 PM  
Printed: 7/19/2006 4:14:00 PM

As performed at the University of Massachusetts, Environmental Engineering Research Laboratory

Prepared by:  
David A. Reckhow      413-545-5392  
18 Marston Hall    [reckhow@ecs.umass.edu](mailto:reckhow@ecs.umass.edu)  
University of Massachusetts  
Amherst, MA 01003

## Table of Contents

Analysis of Total Organic Nitrogen .....	1
Table of Contents .....	1
Standard Operating Procedures .....	3
Analysis of Total Organic Nitrogen .....	3
Scope .....	3
Method Overview .....	5
Detailed Procedures .....	6
Basis for Method .....	6
UMass Detailed Procedures .....	7
Sample Preservation and Pretreatment .....	7
Additional Pretreatment for IN Reduction .....	7
Startup of TOC/TN Analyzer .....	9
Analysis of samples/standards .....	10
TN Analyzer Shut-down procedure .....	12
Data Analysis & QC Reporting .....	13

*UMass Environmental Engineering Program*

Clean Up.....	14
Standard Solutions, Solvents and Supplies.....	15
Preparation of Calibration Standards.....	15
Preparation of Challenge Standards .....	16
Preparation of QC Samples .....	16
Supplies .....	16
Maintenance and Troubleshooting .....	17
Regular Maintenance Activities .....	17
Troubleshooting.....	17
Quality Assurance/Quality Control .....	18
General Approach.....	18
Quality Assurance Objectives.....	19
General Procedures.....	19
Data Quality Indicators.....	19
Precision .....	20
Accuracy .....	20
Method Detection Limit (MDL).....	20
Linearity.....	21
Sampling Custody .....	21
Sample Collection and Storage.....	21
Data Reduction, Validation and Reporting .....	22
Procedures specific to TON Analysis.....	23
General Analytical QC .....	23
Special QC Tests .....	24
Initial Demonstration of Capability (IDC) .....	25
Method Detection Limit (MDL).....	25
QC Criteria .....	25
Literature Cited.....	27
Appendix .....	28
Lee and Westerhoff, 2005   ES&T 39:3:879-884.....	28

## Standard Operating Procedures

# Analysis of Total Organic Nitrogen

This guidance document was prepared to assist research assistants, post-docs and lab technicians in conducting analysis of TON and TN in the UMass Environmental Engineering research laboratories. It aspires to outline our standard operating procedures, as they exist at the present time. It also emphasizes elements of quality control that are necessary to assure high quality data. Please help me keep this document current by alerting me to any errors or long-term changes in methodology or equipment.

Dave Reckhow  
Faculty QC officer for TON analysis

## Scope

This method has been used in the UMass Environmental Engineering Laboratory for analysis of TN and TON in clean water samples (i.e., drinking water, surface waters and uncontaminated groundwaters). It has been found to meet data quality criteria with all waters of these types for which it has been tested. This method should not be used for other media without further validation.

Recoveries using the high temperature oxidation method are generally considered to be complete (i.e., ~100%)<sup>1</sup>. Most tests with purified standards have shown complete recovery within the range of error. The compounds that have been reported to be incompletely recovered are summarized in the table below.

Compounds Reported to be Incompletely Recovered by HT Combustion TN Analysis

Compound	Average Recovery	Source

<sup>1</sup> Recoveries are, of course, all based on an assumed 100% recovery of Potassium Nitrate

Compounds tested at UMass, and known to be Completely Recovered by the Shimadzu TOC-V/TN

Anthranilic acid	4-pyridinealdoxime	Caffeine
3 Aminobenzoic acid	Nicotinic acid	Thymine
4 Aminobenzoic acid	2-hydroxynicotinic acid	Cytosine
2-aminophenol	6-hydroxynicotinic acid	Uracil
Acetanilide	Nicotinamide	3,4,5-Trimethoxybenzamide
1,4-phenylenediamine	Adenine	3,4,5-Trimethoxy phenylacetonitrile
4-hydroxypyridine		Urea

Additional compounds tested by others showing complete recovery<sup>2</sup>

Isonicotinic Acid	Glutaric Acid	Glycine
Glutaric acid	Phenylalanine	Cysteine
Tryptophan		

Other compounds exhibiting complete recovery of N in a Shimadzu TOC-V/TN combustion system<sup>3</sup>

Arginine	Caffeine	EDTA
Histidine	Nicotinic Acid	Sulphathiazole
Thiourea	Thymidine	Urea

Compounds reported by Shimadzu as exhibiting complete recovery in their TOC-V/TN combustion system<sup>4</sup>

Potassium Nitrate	Sodium Nitrate	Ammonium Nitrate
Nicotinic Acid	Urea	

---

<sup>2</sup> From: Wallace et al., 2002

<sup>3</sup> From Alvarez-Salgado et al., 1998

<sup>4</sup> From Walker et al., undated

## Method Overview

Reproduced below is a simple, step-by-step outline of our TON method for quick reference.

**Table 1: Summary of Procedure for TON Analysis**

1. Filter samples if necessary
2. Pretreat samples by Dialysis, if necessary
3. Submit samples for analysis of ammonia (ISE) and nitrate/nitrite (IC)
4. Turn on Analyzer and start log book record
5. Install Autosampler tray
6. Check Gas flow and pressures
7. Check liquid levels
8. Activate furnace from Computer
9. Prepare **calibration standards** (Table 2) and **QC samples** (Table 3)
10. Load samples and standards on tray
11. Record sample identities in log book
12. Start run on computer
13. Verify initial autosampler operation
14. Inspect data and report on QC

**Table 2. Typical Preparation of Calibration Standard**

1. Prepare primary standard containing 100 mg/L of TN in the form of Potassium Nitrate, if needed
2. Prepare intermediate standard (10 mg/L) from the primary standard
3. Prepare working standards (typically 0-2 mg/L) from intermediate standard.

**Table 3. Typical Preparation of QC Samples**

1. Prepare Spiked samples for determination of matrix recovery (laboratory fortified sample matrix). Select 10% of analytical samples and set aside an additional 30 mL aliquot of each.
2. Prepare a continuing calibration check standard at mid-range concentration.
3. Prepare any other QC samples as needed (see Table 6, page 23).

## Detailed Procedures

### Basis for Method

Although there are no standard methods for TON, we use a protocol that is closely aligned with the recommended procedures for operation of the TOC-V/TN as distributed by Shimadzu Corporation. For samples with a high ratio of inorganic nitrogen to organic nitrogen, we use the pretreatment method published by Lee and Westerhoff (2005). Please refer to these documents for all details. However, the analyst should keep in mind that we have occasionally made some specific modifications. Such modifications are itemized below in Table 4.

**Table 4. UMass Protocol Specifics and Departures from Published Methods**

Item #	Step or Apparatus	Published protocol	UMass protocol

## UMass Detailed Procedures

### Sample Preservation and Pretreatment

1. **Filter if Dissolved Organic Nitrogen (DON) is to be measured.**
  - a) Assemble and clean the filtration apparatus
    - May be either glass Millipore-type or in-line syringe filter
  - b) Pre-wash a Whatman GF/C glass fiber filter
    - i. Place filter in housing
    - ii. Rinse with 50 mL deionized water; repeat 5 times
    - iii. If practical, final rinse with 10 mL of sample
  
2. **Drop pH to 2 with one of the following preservatives.<sup>5</sup>**
  - Sufficient  $H_3PO_4$  to bring  $pH \leq 2.0$ 
    - General preservative of choice, unless precipitate forms
  - 100  $\mu L$  of 6 N HCl for every 100 mL of sample.
    - Use this only with pyrolytic system (e.g., Shimadzu TOC-V)<sup>6</sup>
    - Requires that halide scrubber be frequently changed
    - This should be sufficient to lower pH below 2
  
3. **Place aqueous samples in a refrigerator until analysis.**
  - Samples should be analyzed as soon after collection as possible, but under no circumstances should more than 7 days be allowed to elapse.

### Additional Pretreatment for IN Reduction

1. **Determine if Inorganic Nitrogen Reduction is necessary.**
  - a) This requires some prior knowledge of the water samples being tested. When the inorganic nitrogen levels become a substantial background to the total nitrogen (TN), the must be reduced to achieve an acceptable level of accuracy. The total inorganic nitrogen (TIN) is the sum of the nitrate, nitrite and ammonia, as determined by ion chromatography and ion selective electrodes. The default criteria for IN reduction are as follows:
    - $TIN > 0.75 * TN$

---

<sup>5</sup> Nitric acid may be used with UV-persulfate systems, but should be avoided in pyrolytic ones (corrosive  $N_2O_4$  is formed and incompletely removed by scrubber, sulfuric acid should be avoided due to formation of  $SO_3$  in both types of systems (this gas will interfere with  $CO_2$  detection); see Tekmar-Dohrmann application note or Wallace, 2001 [WQTC proceedings])

<sup>6</sup> chloride will scavenge radical oxidants in UV/persulfate systems.

- b) Alternative criteria can be established based on the TON data quality objectives and an error propagation analysis of the TIN and TN data.
- c) If IN reduction is deemed necessary:
  - Set aside an undialyzed sample for measurement of both initial TOC/TN and a UV-Vis spectrum, and
  - Proceed to step 2.
- d) If IN reduction isn't necessary, go directly to IN and TN determination.

**2. Assemble and Load Dialysis System.**

- a) Place 25 mL of each sample into separate 20 mm cellulose ester dialysis tubes
  - Spectra/Por, Spectrum Laboratories, Inc.; nominal MWCO = 100 Daltons
  - Others may be acceptable, especially those with a smaller diameter providing higher surface:volume ratio
  - Tubes must be rinsed with lab grade water for at least 48 hours until the rinse water TOC < 0.1 mg/L
- b) Fill Reservoir with Lab grade water (acceptor solution)
  - Reservoir must be opaque to block out room light, and covered
  - 7L Nalgene tank is acceptable
  - Reservoir must be equipped with a stirrer (either mechanical or magnetic)
- c) Suspend tubes in the dialysis reservoir, and begin stirring

**3. Dialyze for 5 days.**

- a) Run lab fresh grade water through the reservoir for 120 hours (5 days)
  - Use a control device to get the proper flow: either a needle valve attached directly to the DI tap, or a peristaltic pump fed by the DI tap.
  - Flow rate should be set so that it is equal to about half the total sample volume every minute. (i.e., for 9 samples of 25 mL each we would need a flow rate of about 100 mL/min or a total 5-day volume of 720 L)
  - This flow regime has been shown to remove >99% of the nitrate (Lee & Westerhoff, 2005)
- b) At the end of 120 hours, remove the tubes
  - Measure and record a full UV/Vis spectrum of the sample after dialysis for assessment of NOM and nitrate/nitrite loss.
  - Submit the remaining volume for TOC/TN analysis
- c) Drain and rinse the dialysis reservoir



## Startup of TOC/TN Analyzer

### 1. Begin Entry in Log Book

- Record you name, date, number of samples, types of samples
- Also record any problems or unusual performance noted

### a) Turn on the TOC-V/TN Analyzer

- Make sure gas is flowing to the instrument (as observed by bubbling in the IC Reagent Vessel)
- Activate switch on lower left side of analyzer
- Turn on ozone generator if measuring TN
- Wait for indicator to change from “initializing” to “initialized”

### 2. Install autosampler tray.

- Place tray in autosampler, making sure that the positioning slot in the tray is lined up with the positioning pin on the autosampler
- Place turntable cover on top; being sure that the arrow on the cover lines up with the arrow on the autosampler
- Press ASI Initial (F5)
- Wait until autosampler is recognized. The Initial Start screen will reappear when instrument is ready

### 3. Check gas flow and pressure

- Open the knob on the top of the carrier gas cylinder (must be “Ultra Zero” grade air)
- Note first stage pressure on the regulator
  - Cylinder must be replaced when it falls below 500 psi
- Second stage pressure gage on regulator must read 70-85 psi
  - Re-adjust if it doesn't fall in this range
- Open front door of analyzer to reveal gas settings
- Check that the instrument carrier gas pressure gage<sup>7</sup> reads 4-5 kgf/cm<sup>2</sup>
- Check that the instrument carrier gas flowmeter<sup>Error! Bookmark not defined.</sup> reads 150 mL/min
- Verify that the IC reaction vessel is bubbling

### 4. Check Liquid levels

- Verify that the rinse water bottle (located behind autosampler) is full and end of tubing is at the bottom of the bottle
  - If not full, fill with super-Q water and re-adjust tubing
- Verify that the humidifier water level (inside front door of analyzer) is near top white line.
  - If not, unscrew black cap and fill with Super-Q water from a squeeze bottle

---

<sup>7</sup> As long as the pressure from the tank is OK, this setting should be correct. If it is not correct, and the pressure at the tank is within range, there may be some blockages in the gas lines

- c) Verify that dehumidifier drain container (also inside front door of analyzer) is full

**5. Refer to Analysis Section to Turn on the Furnace**

**Analysis of samples/standards<sup>8</sup>**

1. **Prepare standards and QC samples**
2. **Load samples and standard on to autosampler tray**
3. **Place the tray and turntable cover onto the autosampler**
  - make sure the positioning pins and arrows line up
  - See Table 7 (page 24) for typical sample sequence with QC samples
4. **Write down in the log book, sample number and sample identity**
5. **Initiate Sample Run**
  6. Open TOC-V
  7. Select Table Editor (user name TOC)
  8. Open File Tab
    - Select New
    - Select Sample Run
    - Select UMass TN in System Options
    - Save data file template to appropriate location
  - d) Open Insert Tab
    - Select Autogenerate
    - Select appropriate Calibration Curve
    - Next
    - If performing calibration at the beginning of the sample run then select the appropriate option, otherwise Next
    - Choose location of Control Samples
    - Finish
  - e) Sparge/Acid Screen
    - Select OK
  - f) Click on the Lightning Bolt Icon to Connect to the Instrument
    - Use Settings on PC
    - Let instrument initialize
9. **Verify that TOC/V analyzer has warmed up**
  - all status indicators on the “monitor” screen should read “OK”
  - the baseline, shown on the graph, should be flat (although it need not read exactly zero)
10. **Initiate Sample Run**
  - Click on the Stoplight Icon to Begin Sample Run
  - Select Standby Option
  - Confirm Samples/Locations

---

<sup>8</sup> typically requires 10 hours of analyzer time for a run of 30 samples

- Check external acid addition if desired
- START

**11. Verify initial autosampler operation**

- g) After pressing the “START” button, watch as the needles move into the first standard or sample. Verify that the needle is in the right position
- The TOC-V sparges internally

**12. Inspect the first few data points**

- make your first QC report by email (see “Data Analysis” below).

**13. Finish the run and Process peak data when analysis is complete.**

## TN Analyzer Shut-down procedure

1. **If “Finish” was selected in the “ASI Conditions” screen**
  - h) Analyzer will shut off after the last sample was run
  - i) You will see a countdown on the screen
    - “Wait \_\_\_\_ minutes to turn main switch off”
  - j) When countdown is complete, the message “You may now turn main switch off” will appear
  - k) Turn off the analyzer
2. **If “No Change” was selected in the “ASI Conditions” screen.**
  - You can either re-use the analyzer or shut it down manually
  - If you choose to shut down right away, do the following:

## Data Analysis & QC Reporting

- 1) **Data Analysis begins with the first few injections.**
  - a) The analyst must inspect the first few injections to see that:
    - the Super-Q blank is below the threshold for laboratory grade water
    - the first standard exhibits good peak shape
    - reproducibility is acceptable
  
3. **Assessment of Instrument Calibration and Standard Curve**
  - We use least squares best linear fit of the standard peak areas regressed against their known concentrations.
  - Standard curves must also include the zero standard (sometimes called the laboratory reagent blank)
  - Standard curves must be visually inspected for non-linear behavior and the possible presence of outliers
    - When noted, an outlier may be excluded from the calibration curve, after consultation with the graduate QC officer. Removal of an outlier should:
      - Substantially improve the standard curve linearity or correlation
      - Improve agreement with the calibration check standard
      - Bring the regressed slope closer to the expected values based on recent data from the calibration slope control chart
  
4. **Evaluation of calibration check standard and other QC data by the analyst**
  5. This must be done as soon as possible, but no later than 24 hours from the end of the TN run. Compare with quantitative criteria in Table 9.
  - b) Send an email report on the success or failure of these first few injections to the graduate QC officer or his/her designee if he/she is not available. It should include information on:
    - i. Sample types (e.g., field samples from Stamford), field collection date, laboratory treatment date (if any), and analysis date
    - ii. Calibration check performance
    - iii. Blank water (e.g., Super-Q) response
    - iv. Statistics on sample and standard relative standard deviation and number of injections
      - The message must also include the address of the Faculty QC officer in the “cc:” line ([reckhow@ecs.umass.edu](mailto:reckhow@ecs.umass.edu)).
      - The subject line of this email message must simply read “TN QC report”

**6. Validation of QC data**

- a) The graduate QC officer or his/her designee then must compile the analyst's data into the running QC data files, and examine the updated control charts.
- b) The graduate QC officer or his/her designee then must send an email message to the faculty QC officer stating whether the QC data are within control limits, and if they are not, what actions will be taken.
  - Again, the subject line of this email message must simply read "TN QC report".
  - This must be done as soon as possible, but no later than 24 hours from the time of receipt of the detailed QC report (per #3).

**Clean Up**

- All glassware used for TN analysis must be "organic-free"
  - All the vials and any glassware should be rinsed with tap water, let them soak in soapy water overnight and then rinse with RO water and leave them in the acid bath overnight. Then take them out, rinse them with distilled water (3 times) and once with super Q water. Dry them in the oven (The graduated glassware should be placed in the cooler oven)
- Pipette tips may be rinsed with water and super Q water, and let them dry
- Acid baths must be cleaned and refreshed on a weekly basis

## Standard Solutions, Solvents and Supplies

### Preparation of Calibration Standards<sup>9</sup>

#### **1. Prepare Stock Primary Standard.**

- a) Dry about 1 gm of Potassium Nitrate in oven at 103-110°C for 30 min. Cool in desiccator for 20-30 min.
- b) Weigh exactly 0.722 gm using analytical balance. Add to a 1000 mL volumetric flask and fill to mark with Super-Q water.
- c) Result is stock primary standard (100 mg-N/L).
- d) Store in a brown glass bottle. Label with your name, the date and “100 mg-N/L KNO<sub>3</sub> standard.” Note: 1 mL = 0.1 mg TN.
- e) Store in refrigerator. Discard after 2-3 weeks.

#### **2. Prepare intermediate standard.**

- a) Prepare on the day TN samples will be analyzed.
- b) Make a volumetric dilution of the Stock Primary Standard. Pour about 15 mL of the Stock Primary Standard into a beaker. Transfer 10 mL of the stock primary standard with a volumetric pipette to a 100 mL volumetric flask half filled with Super-Q water. Fill to mark with Super-Q water.
- c) Intermediate stock concentration is 10 mg/L TN. Store in refrigerator. Discard after 2 days.

#### **3. Prepare Working standards.**

- a) Prepare 3 working standards that bracket the sample concentrations. For example, for low level TN analysis, a typical calibration curve consists of a 0.5, 0.2 and 0 ppm standard.
- b) Use three 100 mL volumetric flasks. Fill each half way with Super-Q water.
- c) Add the appropriate volume of intermediate stock to each flask (number of mL of intermediate stock = concentration in mg/L of working stock). Fill to mark with Super-Q water. An example for low level TN analysis as shown below.

Working Standard (mg/L)	Volume of intermediate stock added (mL)
10	10
5	5
2	2
0	0

---

<sup>9</sup> typically requires 20 minutes

**Preparation of Challenge Standards**

These are compounds that are relatively resistant to oxidation, especially in low temperature systems.

**1. Prepare Stock Guanine standard.**

**Preparation of QC Samples**

- d) Prepare Spiked samples for determination of matrix recovery (laboratory fortified sample matrix). Select 10% of analytical samples and set aside an additional 30 mL aliquot of each. Add requisite volumes of the intermediate stock so that the TN will approximately double.
- e) Prepare a continuing calibration check standard at a mid-range concentration.
- f) Prepare any other QN samples as needed (see Table 6, page 23)
  - i. Dilute guanine stock to a theoretical TN level of 1 mg/L.

**Supplies**

**Table 5. Summary of Supplies for Organic Nitrogen Analysis**

Item	Catalog #	Approx. Price	Approx # used/run <sup>10</sup>
Pasteur Pipettes	Fisher: 13-678-20A	720/ \$46.10	10
Integrator Paper			
Ink Cartridges	Fisher: 07-684-61	\$15.20 each	
DIUF Water	Fisher: W2-20	\$32.29	Not normally used
H <sub>2</sub> SO <sub>4</sub>	UMass Stockroom		34 mL
Autoamplifier vials			
Autosampler needles			
Potassium Nitrate			
Guanine			

<sup>10</sup> Assuming about 30 samples analyzed



## Maintenance and Troubleshooting

### **Regular Maintenance Activities**

- When measurement sensitivity, repeatability or accuracy has dropped, the catalyst need to be regenerated. This is conducted from the Maintenance screen.
- When regeneration of the catalyst does not correct repeatability and sensitivity problems, the catalyst must be washed or replaced.
- Replace the Ultra Zero air when tank pressure falls below 500 psi.
- Replace the CO<sub>2</sub> absorber annually.
- Replace the halogen adsorber before the colored area reaches the front end.
- Replace the rinse water pump head after 300 hours of use.

### **Troubleshooting**

Past difficulties have resulted from:

- Improper positioning of the ASI unit with respect to the analyzer.
  - If the ASI is pushed too far back, it cannot properly collect samples, nor is it obvious that this is the source of the problem

For further assistance:

- Consult Shimadzu manuals
- Call Shimadzu technicians (800 number)

## Quality Assurance/Quality Control

### General Approach

Quality assurance is an essential and integral part of a research study. The purpose of any QA plan is to insure that valid and reliable procedures are used in collecting and processing research data. The procedures outlined are designed to eliminate or reduce errors in experiments, sample preparation and handling, and analytical methods. Emphasis must be given throughout one's lab work to incorporate the plan into the research project by all research personnel.

Any equipment and experimental procedures that are used to provide numerical data must be calibrated to the accuracy requirements for its use. Records are to be kept of all calibrations. Calibration schedules are generally established for all aspects of physical and chemical measurements and these must be strictly followed. Physical standards and measuring devices must have currently valid calibrations, traceable to national standards. Most chemical standards are acquired from commercial suppliers, and they should be of the highest purity available. When necessary, unavailable standards should be synthesized using the best methodology available.

As a general rule, experiments should be replicated to assure reproducibility. All data reported should include a statement of its uncertainty, and the means for the determination and assignment of such limits. Standard reference materials are used for this purpose where possible. Statistically established confidence limits and an analysis of sources of systematic error are to be used in the absence of experimental demonstration of limits of inaccuracy.

All data will be subject to review by the faculty QC officer before release. The analysts involved will certify reports as well as all who review them. All analysts and QC officers must attest that the data and associated information contained in the report are believed to be correct and that all quality assurance requirements have been fulfilled, unless exceptions are approved and noted. Careful and detailed laboratory records will be maintained by each analyst, including source of reagents, meticulously detailed procedures (referring to an SOP, and any departures or clarifications), instrument and conditions of analysis, failed experiments, etc. Data output will be archived.

Regular meetings will be held to review the results and project progress, and to plan further experiments. The results will be analyzed promptly and summarized by means of internal reports or formal reports for external review. The experimental and analytical procedures will be reviewed for their performances and changes will be made as necessary. The quality assurance program as described in this document must be strictly observed.

## **Quality Assurance Objectives**

The precision, accuracy and method detection limits will be evaluated, and where there are existing methods, held within the control limits set forth in the accepted references (e.g., APHA et al., 1999; USEPA-EMSL, 1990; ASTM, 1994). In addition to the analysis of sample replicates, a minimum of 10 percent of the time is typically involved in analytical determinations that are devoted to quality control. The precision of each test is determined through analysis of sample replicates. These are commonly presented in the form of control charts (e.g. Section 1020B of APHA et al., 1999).

The accuracy of each analysis will be determined by measuring spike recoveries in the matrix of interest. The relative errors will be calculated and will be considered acceptable if they fall within the control limits determined for the particular test. For new methods developed at UMass or for modifications of existing methods, we will have to establish criteria on acceptable control limits. In general, a test will not be deemed useful if its precision or accuracy is found to be equal to or greater than 20% of the highest values observed. Where possible, external performance standards will also be run. This serves as a measure of accuracy both for the analysis and for standard preparation.

Data generated by the QA program will be incorporated into a Quality Control (QC) archive that is used to monitor the fluctuations in precision and accuracy so that chance or assignable causes of error can be determined. Control charts such as X-charts for simple successive samples or cumulative sum techniques may be employed to record both precision and accuracy data (Taylor, 1987).

## **General Procedures**

General sample collection and handling will be in accordance with the guidelines of Section 1060 of Standard Methods (APHA et al., 1999). All previously established analytical methods used in laboratory research will follow approved methods in the standard compilations (e.g., APHA et al., 1999; USEPA-EMSL, 1990, or ASTM, 1994).

Reagent grade chemicals or higher quality when needed will be used throughout the research. Super-Q water (purified by reverse osmosis, deionization, and carbon adsorption) will be used for preparation of reagents, sample blanks, and dilution water. Where necessary, this water will be further purified using batch UV irradiation. All glassware used in the experiments and in analytical analyses will be thoroughly cleaned with a chromium-free sequence of detergent, oxidant and acid to prevent interferences from trace contaminants.

## **Data Quality Indicators**

A wide range of data quality indicators are normally calculated for the purpose of assessing method performance. Some of these are defined below.

### **Precision**

Precision may be expressed as the relative percent difference (RPD) from duplicate measurements ( $C_1$  and  $C_2$ ) of the same sample:

$$RPD = \frac{|C_1 - C_2| \times 100\%}{(C_1 + C_2)/2}$$

When three or more replicates are available, the relative standard deviation (RSD) should be used:

$$RSD = \left( \frac{s}{\bar{y}} \right) \times 100\%$$

where the standard deviation (s) is determined from:

$$s = \sqrt{\frac{\sum_{i=1}^n (y_i - \bar{y})^2}{n-1}}$$

### **Accuracy**

Accuracy is best assessed by analysis of a standard reference material (SRM) prepared in the matrix of interest. It is quite rare that such materials are available, so two possible compromises may be used instead. These are the laboratory-prepared matrix spikes, and the independent SRM prepared in a standard matrix. One or both may be analyzed and the percent recovery (%R) calculated as a measure of accuracy.

$$\%R = \left( \frac{S - U}{C_{sa}} \right) \times 100\%$$

where:

S = measured concentration in spiked aliquot

U = measured concentration in unspiked aliquot

$C_{sa}$  = actual concentration of spike added

$$\%R = \left( \frac{C_m}{C_{sm}} \right) \times 100\%$$

and:

$C_m$  = measured concentration of SRM

$C_{sm}$  = actual concentration of SRM

### **Method Detection Limit (MDL)**

The method detection limit (MDL) will be defined as:

$$MDL = s_7 \bullet t_{(n-1, 1-\alpha=0.99)}$$

where:

$s_7$  = standard deviation of 7 replicate analyses where the mean is no more than 10 times the MDL

$t_{(n-1, 1-\alpha=0.99)}$  = Students' t-value for a one-sided 99% confidence level and a standard deviation estimate with n-1 degrees of freedom.

### **Linearity**

The calibration curve linearity (L) is defined as the ratio of the slope using the highest standards ( $S_U$ ) divided by the slope determined from the lowest standards ( $S_L$ ) as follows:

$$L = \frac{S_U}{S_L}$$

The highest standards shall be all those that fall within the top 50% of the calibration range including the 50% standard if it is used. If only one standard falls within that range, the  $S_U$  shall be calculated based on the top two standards. The lowest standards are all those that fall within the bottom 50% of the calibration range including the 50% standard if it is used. Least squares linear regression is used to determine slopes.

### **Sampling Custody**

In most cases analyses will be performed immediately upon return from the field or after preparation of samples in the laboratory. Problems with sample custody are minimized, because the same people who receive (or sometimes, collect) the samples also analyze them. In general sample collection, handling, and preservation will be in accordance with the guidelines of Section 1060 of Standard Methods (APHA et al., 1999). All samples must be fully labeled with the sample identity, date, and name of researcher.

### **Sample Collection and Storage**

Samples are collected and stored in clean borosilicate (e.g., Pyrex, Kimax) glass containers. Containers must be capped with either Teflon-lined septa or ground glass stoppers. Exceptions are made for large volume samples which may be stored and shipped in clean polyethylene carboys. Glass containers are cleaned with detergent, followed by 5% sulfuric acid soak, and final rinsing with reagent-grade water. These containers are dried in a 150°C oven.

Samples for TOC analysis must be acidified and kept in a refrigerator from the time of collection to the start of analysis. Some organics are biodegradable, so care must be taken to minimize this type of loss.

## **Data Reduction, Validation and Reporting**

To ensure the accuracy and permanency of collected data, all research data are recorded with permanent ink in bound notebooks and all QC data (precision, accuracy) are recorded in instrument log notebooks. Summary QC graphs and tables are reviewed at least quarterly by the Faculty QC officer to observe noteworthy trends or inconsistencies. These are maintained in loose leaf notebooks for subsequent use in preparing progress reports, final reports, and theses. Major concerns and conclusions are reported to the external Project Officer via the progress reports.

Pages from the laboratory data books are regularly duplicated so that a file copy of raw data can be placed in safe storage in the event that the book is lost or destroyed. At the end of the project, all bound data books and any loose leaf data will be stored by the project team for at least ten years. Summary data files will be put on magnetic media so that statistical analysis of the data can be done. Our laboratory has several personal computers that can be used for this purpose.

## Procedures specific to TON Analysis

### General Analytical QC

Many types of QC procedures are required as indicated under standard method 5310B for TOC, and most of these have been adapted to TON analysis. The guidelines below are prepared assuming that samples are run in groups, whereby a “daily” frequency refers to once every day that the analytical method is being used.

**Table 6. Summary of QC Elements as Applied to TON Analysis**

Types of Samples or Standards	§ in Std. Meth.	Purpose	Frequency	Timing	QC data
Laboratory Performance Check Standard (LPC)		To establish basic analyzer performance	1 standard per day	Beginning of each day	Noise, response
Initial Demonstration of Capability (IDC)		To show that an analyst’s technique and equipment are adequate for TON analysis	4-7 standards when first learning method, otherwise not done		Mean % recovery and standard deviation
Method Detection Limit (MDL)		To determine the lowest concentration level that the analyst can report	7 standards run		MDL and EDL
Laboratory Reagent Blank (LRB)		Test lab conditions and quench for interferences	1 per day	Beginning of day	Response
Field Reagent Blank (FRB)		Test all field conditions for interferences	1 per day, if sampling occurred outside of the lab	mid day	Response
Spiked sample, or Laboratory Fortified Sample Matrix (LFM)		To test recovery in the sample matrix	1 for every 10 samples <sup>11</sup>	Mixed throughout day	% recovery, mean and standard deviation
Calibration Standards		To provide a basis for determining the concentrations in unknowns	3 levels including zero		Calibration curves including slopes and intercepts
Continuing Calibration Check		To verify the accuracy of the calibration	Usually one per day	mid-day	

<sup>11</sup> These are generally done using the calibration compound (e.g., KHP). Spike recovery studies may also be done with a “challenge” compound.

Standards (CCC) <sup>12</sup>		standards			
Challenge standard <sup>13</sup>		To test the method with a “difficult” compound	Irregular: e.g., when catalyst is changed		% recovery
Unknowns or “samples”		This is what you really want to measure	As many as desired	Mixed throughout day	Relative standard deviations

Table 7 shows a recommended sequence for a typical run of about 16 samples. The first two samples require immediate attention, as they are simple indicators of unacceptable QC. When these detector response traces show abnormally elevated noise levels or drifting baseline, the operator must intervene before proceeding. The problem must be diagnosed, solved and the sequence restarted at vial #1.

**Table 7: Typical Vial Sequence for TN Analysis**

Vial #	Sample type	QC objectives
1	Water Blank	To check on IN condition
2	Zero Standard	To check for gross contamination of water or lab environment
3-5	Remaining calibration standards	Calibration
6-13	Analytical Samples	
14	Spiked sample	Spike recovery
15	Calibration check standard	Check on accuracy of calibration samples
16	Calibration standard (one intermediate)	Final check to verify that calibration hasn't changed during run

### Special QC Tests

The following are “special” tests that are not part of the normal QC protocol. They are used when first learning this analytical method (e.g., IDC), and they may be used when there are suspected problems or there is a need for method performance evaluation.

<sup>12</sup> Prepared from the calibration solution (eg., KNO<sub>3</sub> stock)

<sup>13</sup> these are compound that are not as easily oxidized, such as caffeine, and guanine. Sometimes standard solutions of humic substances can be used.



### **Initial Demonstration of Capability (IDC)**

This is normally performed by each analyst when he/she is first learning to measure TON.

#### ***Procedure***

- a) Prepare a full standard curve
- b) Perform an MDL test for TON
- c) Prepare 3 sets of samples with matrix spikes, and evaluate spike recovery in 3 different matrices

#### ***Performance***

IDC results are expected to meet the LFM and MDL QC criteria from Table 8.

### **Method Detection Limit (MDL)**

#### ***Procedure***

- a) Secure a 1000 mL sample with a low level of TON (preferably below 0.3 mg/L). This may be a bulk sample of tap water diluted 50/50 with Super-Q water.
- b) Separate this into 7 aliquots of 30 mL each.
- c) Analyze each for TON on the same day.
- d) Determine MDL based on the standard deviation of these 7 sets of measurements (refer to section on: MDL calculations on page 25 of this document)

#### ***Performance***

We expect that MDLs for TON will be below 0.1 mg/L. Shimadzu cites 0.01 mg/L as an MDL in their laboratory under ideal conditions.

### **QC Criteria**

Quality control data must be analyzed as soon as possible. The best practice is to have the QC data tabulated and evaluated as the run is underway. However, it is recognized that there will be times when this is impossible. QC and calibration data must always be analyzed and reported within 24 hours of completion of a run (see section on

Data Analysis & QC Reporting, page 13). Quantitative criteria (Table 8) must be applied, and violations must be immediately reported to the faculty QC officer. The graduate and faculty QC officer along with the analyst will then work out a plan for returning the analysis to acceptable levels of QC. Table 8 lists some typical corrective action, however the actions taken may differ depending on the particular circumstances. Excursions from QC criteria can be quite complex, and many analytical characteristics and conditions must be considered before a decision can be made on the most effective steps to be taken.

In several cases, quantitative criteria are based on long term trends, and these must be monitored by means of appropriate control charts. Calibration check standards, spike recovery and blank water values are documented over time in this way. All summarized QC data (tabular and graphical) must be kept in a notebook in room 308 of Elab II. A duplicate set must be deposited with the faculty QC officer (D. Reckhow).

**Table 8: Quantitative Criteria for Judging Data Acceptability**

Types of Samples or Standards	Frequency	Timing	QC data Acceptance Criteria	Typical Corrective Action
Spiked sample, or Laboratory Fortified Sample Matrix (LFM)	1 for every 10 samples	Mixed throughout day	Mean % recovery = 80%-120%	<ul style="list-style-type: none"> <li>❖ Re-run matrix spikes</li> <li>❖ Re-examine entire run for errors</li> <li>❖ Possibly change SOP</li> </ul>
Calibration Standards <sup>14</sup>	4 levels including zero		Calibration slopes = $\pm 30\%$ of long-term average	<ul style="list-style-type: none"> <li>❖ Run new set of standards</li> <li>❖ Prepare new anion stock</li> <li>❖ Examine TOC analyzer for problems, needed maintenance</li> </ul>
			Calibration linearity must not fall below 0.5	❖
Continuing Calibration Check Standards (CCC) <sup>15</sup>			Calculated Conc. = $\pm 20\%$ of expected value	<ul style="list-style-type: none"> <li>❖ 1. Prepare new calibration check standard</li> <li>❖ 2. Prepare new standard curve based on new stock</li> </ul>
Method Detection Limit (MDL)	At IDC and once per year		Below 0.1 mg/L	<ul style="list-style-type: none"> <li>❖ Examine sample handling</li> <li>❖ Check instrument of baseline noise level</li> </ul>
			MDL for DOC should be no more than 0.02 mg/L above MDL for TOC	<ul style="list-style-type: none"> <li>❖ Re-evaluate filtration procedures</li> </ul>
Lab Reagent Blank	Daily		$\leq 0.02$ mg/L	<ul style="list-style-type: none"> <li>❖ Seek alternative reagent water</li> <li>❖ Change cartridges</li> </ul>

<sup>14</sup> Prepared from the currently-used calibration stock (less than 1 month old)

<sup>15</sup> Prepared from the previously-used calibration stock II


Unknowns or "samples"	As many as desired	Mixed throughout day	RSD or RFD for replicate analyses $\leq 20\%$	❖ Re-run samples and/or discard outliers <sup>16</sup> until precision can be brought under control
			Estimated concentration in unknowns must not exceed highest standard	<ul style="list-style-type: none"> <li>❖ Re-run samples with higher level standards</li> <li>❖ Dilute and re-run samples</li> <li>❖ If within 150% of max standard, concentrations may be flagged as tentative</li> </ul>

## Literature Cited

- Alvarez-Salgado, Xose, and Axel E.J. Miller, (1998) *Marine Chemistry* 62: 325-333.
- Egeberg, P.K., Eikenes, M., and Gjessing, E.T. (1999) Organic Nitrogen Distribution in Nom Size Classes. *Environment International* 25 (2-3), 225-236.
- Hagedorn, F., Schleppe, P., Waldner, P., and Fluhler, H. (2000) Export of Dissolved Organic Carbon and Nitrogen From Gleysol Dominated Catchments - the Significance of Water Flow Paths. *Biogeochemistry* 50 (2), 137-161.
- Lee, W.T. and Westerhoff, P. (2005) Dissolved Organic Nitrogen Measurement Using Dialysis Pretreatment. *Environmental Science & Technology* 39 (3), 879-884.
- Sowden, F. J., Griffith, S. M., and Schnitzer, Morris. The Distribution of Nitrogen in Some Highly Organic Tropical Volcanic Soils. *Soil. Biol. Biochem.* 8, 55-60. 1975.
- Spyres, Georgina, Malcolm Nimmo, Paul J. Worsfold, Eric P. Achterberg and Axel E. J. Miller (2000) Trends in Analytical Chemistry 19(8)498
- Walker, Karmel R., Lisa Stojowski, and Robert H. Clifford, (undated) Total Nitrogen Analysis: A New Perspective on TOC, Shimadzu website ([www.shimadzu.com](http://www.shimadzu.com))
- Wei, Ying, Jim Reavis, Curtis Meadows, Lillian Woung and Jim Greenlee (1998) WQTC Proceedings
- Westerhoff, P. and Mash, H. (2002) Dissolved Organic Nitrogen in Drinking Water Supplies: a Review. *Journal of Water Supply Research and Technology-Aqua* 51 (8), 415-448.

<sup>16</sup> using Dixon's Q Test, or some logic test (e.g., monotonic increase with timed data series; unlikely sudden changes in natural system).

# Appendix



Lee and Westerhoff, 2005

ES&T 39:3:879-884