

Waters

Eric S. Grumbach

Thomas E. Wheat

Jeffrey R. Mazzeo


Diane M. Diehl

Increased Throughput and Sensitivity Obtained
with UPLC™ Columns Packed with
1.7 μm Particles for High Resolution, Ultra Fast
Stability-Indicating Method Development





For Complete ☼ Confidence

- Pharmaceutical scientists must develop assay procedures that completely identify and measure all degradation products of an active pharmaceutical ingredient.
- Due to its ability to separate degradation products, excipients and process impurities from active ingredients, HPLC has become the analytical tool of choice for stability-indicating assays.
- However, there is always the requirement to achieve better resolution to ensure complete characterization of the degradants.
- At the same time, improvements in sensitivity to detect trace level components and improved sample throughput need to be addressed.
- These assays can benefit from utilizing sub-2 μm particulate columns to improve resolution for critical pairs or maintain existing resolution while improving sample throughput.
- In this study, we examine this approach, applying Ultra Performance LC™ (UPLC™) to the degradants of the antifungal terbinafine.



Faster, More Sensitive Methods

**Faster, More Sensitive,
Higher Resolution Methods**



Choices available:

- Increased **speed** and sensitivity with the same resolution
- Increased **sensitivity** and speed with resolution
- Increased **resolution** and sensitivity at the same speed




- A new class of separation science
 - Based on chromatography columns with very small particles
 - Based on instruments designed to take advantage of the small particles
- Provides improved Resolution, Speed, and Sensitivity
- Suitable for chromatographic applications in general
 - Appropriate for improving existing methods
 - Appropriate for developing new methods

- Sub-2 µm particles
 - Porous for optimum loadability
 - Bridged hybrid particle* required for both high strength and outstanding chromatographic performance
 - Innovative sizing technology for narrow particle size distribution
- Column hardware
 - New frit technology* to retain particles
 - Fittings optimized for high pressure operation
- Packing technology
 - New column packing processes to optimize stability
- eCord™
 - New information chip to store column history


*Waters Patented Technology
No. 6,686,035 B2

Introducing 2nd Generation Hybrid: Bridged EthylSiloxane/Silica Hybrid Particles



1.7 μm BEH particle provides:

- Improved Strength
- Improved Efficiencies
- Improved Peak Shape
- Wider pH Range



Available Ligands:

- ACQUITY UPLC™ BEH C₁₈
 - Straight chain alkyl C₁₈
- ACQUITY UPLC™ BEH Shield RP₁₈
 - Embedded polar group (carbamate)
- ACQUITY UPLC™ BEH C₈
 - Straight chain alkyl C₈
- ACQUITY UPLC™ BEH Phenyl


Why Multiple Ligands?:

- Changes in hydrophobicity
- Changes in silanol activity
- Changes in hydrolytic stability
- Changes in ligand density
- Changes in selectivity

- Specialized instrumentation was designed to meet the requirements of 1.7 μm packings.
- These requirements include:
 - low system volumes (<140 μL to minimize dispersion),
 - high pressure fluidic modules (up to 15,000 PSI),
 - high speed optical detectors (capable of 40 Hz) and
 - mass detectors (5000 Da/sec).
- The integration of both chemistry and instrumentation leads to the development of ultra-fast, sensitive, high resolution chromatographic methods.

Sample Preparation:

- Terbinafine HCl was forcefully degraded with 8.0 N hydrochloric acid
- A 10 mg/mL solution was stirred in a 60 °C water bath for 60 minutes
- A 1 mL aliquot was then neutralized with sodium hydroxide.
- This solution was then diluted with 10 mL of water and 15 mL of acetonitrile for analysis on the UPLC™ system.



Terbinafine
mw = 291.43

Instrument:

ACQUITY UPLC™ with TUV optical detector

Chromatographic Conditions:

Column: C₁₈ 2.1 x 150 mm, 5.0 µm

Mobile Phase A: 20 mM ammonium bicarbonate pH 10.0

Mobile Phase B: acetonitrile

Flow Rate: 0.2 mL/min

Isocratic: 73% B

Injection Volume: 4.0 µL

Sample Diluent: 50 ACN with 10 mM NH₄HCO₃ pH 10


Temperature: 30 °C

TUV Optical Detector Settings:

Wavelength: 214 nm

Sampling Rate: 40 Hz


Time Constant: 0.1 seconds



Choices available:

- Increased speed and sensitivity with the same resolution
- Increased sensitivity and speed with resolution
- Increased resolution and sensitivity at the same speed





L/dp

RATIO

$$\frac{150\text{mm}}{5\mu\text{m}}$$

=

30,000

$$\frac{50\text{mm}}{1.7\mu\text{m}}$$



=

29,500

$$\frac{100\text{mm}}{1.7\mu\text{m}}$$

=

58,820



(F ↑ 2X, L ↓ 3X)

Speed ↑ 6X

Sensitivity ↑ 2.5X

Rs = 1X

- Greater sensitivity
 - Inject same volume onto a 1.0 mm i.d. column packed with 1.7 μm particles
- Greater resolution
 - Use a 100 mm column packed with 1.7 μm particles
($L/d_p = 58,820$)

- Greater sensitivity

- Inject same volume onto a 1.0 mm i.d. column packed with 1.7 μm particles


- Greater resolution

- Use a 100 mm column packed with 1.7 μm particles
($L/dp = 58,820$)

Choices available:

- Increased speed and sensitivity with the same resolution
- **Increased sensitivity and speed with resolution**
- Increased resolution and sensitivity at the same speed






Speed \uparrow 6X

Sensitivity \uparrow 8X



$R_s = \downarrow$ 1.55X



•Same amount
 injected on
 column results in
 loss in R_s

•This is due to
 column overload

Improved Throughput and Sensitivity: 1.0 mm i.d. Column Scaled Injection Volume



**Scaled injection
volume**

4.0 µL → 0.9 µL

- Greater sensitivity

- Inject same volume onto a 1.0 mm i.d. column packed with 1.7 μm particles

- Greater resolution

- Use a 100 mm column packed with 1.7 μm particles
($L/d_p = 58,820$)

L/dp

RATIO

$$\frac{150\text{mm}}{5\mu\text{m}} = 30,000$$



$$\frac{50\text{mm}}{1.7\mu\text{m}} = 29,500$$

$$\frac{100\text{mm}}{1.7\mu\text{m}} = 58,820$$

Choices available:

- Increased speed and sensitivity with the same resolution
- Increased sensitivity and speed with resolution
- **Increased resolution and sensitivity at the same speed**





Speed = 1X

Sensitivity \uparrow 2X

Rs = \uparrow 1.8X

Benefits of applying UPLC™ to stability-indicating assays:

- Develop **faster** and more sensitive assays while maintaining existing resolution (2.1 mm i.d.)
 - **Speed** ↑ **6X**
 - **Sensitivity** ↑ **2.5X**
- Develop highly **sensitive** and faster assays (1.0 mm i.d.)
 - **Sensitivity** ↑ **8X**
 - **Speed** ↑ **6X**
- Develop sensitive, higher **resolution** assays within the same analysis time
 - **R_s** = ↑ **1.8X**
 - **Sensitivity** ↑ **2X**