



ACQUITY UPLC[™] in the Chromatography Lab: Analytical Considerations

Ultra Performance LC[™] delivers improved speed, sensitivity and resolution that can benefit many liquid chromatography applications in today's laboratory. We will discuss analytical performance along with considerations and recommendations for successful implementation of UPLC[™] in the chromatography laboratory.

©2005 Waters Corporation

Why Advance the Science and Technology of LC?

Challenges in the LC Lab

- Chromatography
 - Resolution
 - Reproducibility
 - Ruggedness
 - Sensitivity
 - Speed
 - Detector optimization

- Business
 - Complex samples
 - Regulatory considerations
 - Increased # samples
 - Time pressures
 - Cost savings

Can you meet today's challenges?

HPLC: 7 Analgesics in 25 minutes

©2005 Waters Corporation

Gradient Table

Time	Flow	%A	%В	%C	%D	Curve
	1.00	98.0	2.0	0.0	0.0	
25.00	1.00	75.0	25.0	0.0	0.0	6
25.01	1.00	98.0	2.0	0.0	0.0	6
60.00	0.05	98.0	2.0	0.0	0.0	6

Components

Peak #1 – Acetaminophen

Peak #2 - 2-Acetomidophenol

- Peak #3 Caffeine
- Peak #4 Acetanilide
- Peak #5 Acetylsalicylic Acid
- Peak #6 Salicylic Acid
- Peak #7 Phenacetin



UPLC[™]:7 analgesics in 25 seconds

©2005 Waters Corporation



Minutes

Why Change from HPLC to UPLC™? Time and solvent savings

©2005	Waters Corporation					-		
			Column Volume (X5)	System Volume (x3)	Re- equilibr ation time (min.)	Total cycle time (min.)	# of injections in 24hours	Solvent consumed per injection
	HPLC	1.0	1.66mL	620uL	10.2.	35.2	41	35.2mL
	4.6 x 100 5μ 25 min. run time	mL/min.	(8.3mL)	(1.86mL)				
	UPLC™	1.3	0.20mL	120uL	1.04	1.45	993	1.9mL
	2.1 x 50 1.7μ	mL/min	(1.0mL)	(0.360mL)				
	0.41 min run time							

Total Liters of solvent for 1000 injections: HPLC: 35.2L Vs. UPLC[™]: 1.9L

Why Change from HPLC to UPLC™? More information about your sample





Why Change from HPLC to UPLC™? Improved sensitivity

©2005 Waters Corporation

100 pg/mL 1.72 Mass1 > Mass2 2.08 796 Rt = 2.08 min, Customer Compound A, 100pg/mL, 10µL injection. **HPLC** Phenomenex column 2.0x30mm, 4um Peak Width = 6 secs 8 Time 1.00 2.00 3.00 100 pg/mL 0.83 Mass1 > Mass2 Rt = 0.83 min, Customer Compound A, 1.15e4 100pg/mL, 10µL injection. 0.54 ACQUITY UPLC[™] Column C18 2.1x50mm, 1.7µm Peak Width = 1.8 secs **UPLC** ~ Time 0.30 0.60 1.00

Method Comparison 100pg/mL, 10µL injection

Speed, Sensitivity, Resolution

©2005 Waters Corporation



©2005 Waters Corporation

HPLC Today The state of the science

 The quality of the separation in HPLC is both driven and limited by column chemistry

 Optimized chemistry can bring speed, sensitivity and resolution benefits to the chromatographic process

 Particle size and selectivity are primary contributors to column performance

 Today's classical HPLC instruments cannot maximize the reality of "state of art" smaller particles

©2005 Waters Corporation



van Deemter plot for classical HPLC

Smaller Particles The enabler of productivity

- Smaller particles provide:
 - increased efficiency
 - maintain efficiency over a wider linear velocity
 - ability for both added resolution and increased speed of separation
 - Particles are central to the quality of the separation

Fundamental Resolution Equation At constant column length

Retentivity

©2005 Waters Corporation

$$\mathsf{Rs} = \frac{\sqrt{\mathsf{N}}}{4} \left(\frac{\alpha - 1}{\alpha} \right) \left(\frac{\mathsf{k}}{\mathsf{k} + 1} \right)$$

System Selectivity Efficiency

•In UPLC[™] systems, N (efficiency) is the primary driver

- •Assuming the Selectivity and retentivity are the same as in HPLC
- •Then, Resolution, Rs, is proportional to the square root of N

$$Rs \propto \sqrt{N}$$

And, Efficiency (N), is inversely proportional to Particle Size , dp

So:
$$N \propto \frac{1}{dp}$$

Therefore if: dp \downarrow 3X, Then: N \uparrow 3X, Rs \uparrow 1.7X

Smaller Particles The enabler of productivity



Summary of Key System Attributes (at constant column length)

©2005 Waters Corporation

	Resolution Improvement	Speed Improvement	Sensitivity Improvement	Back Pressure
1.7 vs. 5 µm particles	1.7X	3X	1.7X	27X
1.7 vs. 3 µm particles	1.3X	2X	1.3X	6X

Influence of Particle Size on Chromatography



Isocratic methods with 2.1 x 50 mm columns

Optimal Linear Velocity on ACQUITY UPLC[™]: Isocratic Caffeine Metabolites Separation

©2005 Waters Corporation

HPLC Non-Optimal Linear Velocity





ACQUITY UPLC[™] Technology Ultra performance by design



©2005 Waters Corporation

How do we take advantage of smaller particles?

- Pressure tolerant particle/column
- High pressure fluidic modules
- Holistic design, low system volume, optimized flow paths
- Reduced cycle time autosampler with minimum carryover
- High speed detectors; optical and mass
- Controlled and coordinated system interaction
 - Software designed for system integration
 - Comprehensive diagnostic suite

©2005 Waters Corporation

ACQUITY UPLC[™] Chemistry Proposition

Advancing the Science

- Requires innovation in every aspect of column:
 - New patented bridged hybrid particle required for pressure tolerance and outstanding chromatographic performance
 - Sub 2µm particles
 - Porous for optimum mass transfer
 - Column hardware
 - New patented frit technology to retain particles
 - New end fittings for high pressure/low dispersion operation
 - Packing technology
 - New column packing processes to optimize stability
 - eCord[™]
 - New information chip to store column history

Reversed-Phase Column Selectivity Chart

©2005 Waters Corporation 3.6 Waters Spherisorb® S5 P 3.3 **ACQUITY UPLC™** 3 (In [lpha] amitriptyline/acenaphthene) **BEH C18 ACQUITY UPLC™** 2.7 Waters Spherisorb® S5CN **BEH Phenvl** Nova-Pak® CN HP Waters Spherisorb®/ODS1 2.4 Inertsil® Ph-3 Resolve® C18 Hypersil® Phenyl 2.1 Waters Spherisorb® ODS2 µBondapak™ C18 1.8 YMC-Pack[™] YMC J'sphere[™] Hypersil® CPS Cyano Inertsil® CN-3 1.5 Nucleosil® C18 Phenvl ODS-L80 Nova-Pak® Phenyl YMC J'sphere™ ODS–M80 1.2 Hypersil® BDS Phenyl 🌰 YMCbasic™ Chromolith[™] Nova-Pak® YMC J'sphere™ ODS–H80 XTerra ® **RP-18** C18 0.9 Phenvl Luna ™**♦** YMC-Pack[™] ODS–AQ[™] Nova-Pak® YMC-Pack[™] CN YMC-Pack[™] Pro C4[™] Phenvl Hexvl Atlantis[®] dC18 YMC-Pack™ YMC-Pack™ ODS-A™ 0.6 ACT Ace® C18 Symmetry® C8 **ACQUITY UPLC™** Zorbax® XDB C18 XTerra[®] MS C8 Inertsil® ODS-3 **BEH C8** 0.3 Pro C18™ ♦ SunFire ™ C18 Luna ® SunFire [™] C8[°] XTerra[®] MS C18 C18 Symmetry® C18 SymmetryShield™ RP8 (2)0 Zorbáx® SB C18 XTerra® RP18 ♦ SymmetryShield™ RP18 -0.3 XTerra® RP8 YMC-Pack[™] PolymerC18[™] -0.6 -1.5 -0.5 0.5 1.5 2.5 3.5 (In [k] acenaphthene) **ACQUITY UPLC™** 012005 Shield RP18

ACQUITY UPLC[™] Column Selectivity Comparison for Analgesics

©2005 Waters Corporation



- 1. Acetaminophen
- 2. 2-Acetomidophenol
- 3. Caffeine
- 4. Acetanilide
- 5. Acetylsalicylic Acid
- 6. Salicylic Acid
- 7. Phenacetin

ACQUITY UPLC[™] BEH C18 1.7µm ACQUITY UPLC[™] BEH RP18 1.7µm ACQUITY UPLC[™] BEH C8 1.7µm ACQUITY UPLC[™] BEH Phenyl 1.7µm

©2005 Waters Corporation

Column Lifetime: >8,500 Injections



Column Lifetime Customer Example

©2005

.

Isocratic aging experiment

- Cycle: low pressure runs/injections (data shown) → 200/400/600 etc. high pressure runs/injections (13,000 psi) → low pressure runs/injections
- This aging cycle was repeated
- Columns lasted over 4,000 high pressure injections

Conditions (Isocratic Slow Flow/Low Pressure)

Column:	ACQUITY UPLC [™] BEH C18 (2.1 x 50 mm, 1.7 μm)
	ACQUITY UPLC [™] BEH Shield RP18 BEH (2.1 x 50 mm, 1.7 μm)
Part Number:	186002350
Mobile Phase A:	0.1% HCOOH in Waters
Mobile Phase B:	0.1% HCOOH in ACN
Flow Rate:	0.5 mL/min (~6,000 psi)
Condition:	67% A / 33% B
Injection Volume:	5.0 μL
Sample Concentration:	208 μ g/mL ketorolac tromethamine and 13 μ g/mL naproxen in MeOH
Temperature: 65°C	
Detection: UV 225 nm (Sample	e Rate: 10 Hz, Filter: Normal)

Conditions (Isocratic Fast Flow/High Pressure)

Columns:	ACQUITY UPLC [™] BEH C18 (2.1 x 50 mm, 1.7 μm)		
	ACQUITY UPLC [™] BEH Shield RP18 (2.1 x 50 mm, 1.7 μm)		
Mobile Phase A:	0.1% HCOOH in Water		
Mobile Phase B:	0.1% HCOOH in ACN		
Flow Rate:	1.250 mL/min (Pressure ~13,000 psi)		
Condition:	62% A / 38% B		
Injection Volume:	5.0 μL		
Sample Concentration:	208 $\mu\text{g}/\text{mL}$ ketorolac tromethamine and 13 $\mu\text{g}/\text{mL}$ naproxen in MeOH		
Temperature: 65°C			
Detection: UV 225 nm (Sample	e Rate: 10 Hz, Filter: Normal)		

Column Lifetime Customer Example



Chromatograms obtained under low (~6000 psi) isocratic conditions



©2005 Waters Corporation

- What about chemical stability of the column at the extremes of pH, temperature and pressure?
- A 2 x 2 experimental design (at UPLC[™] pressures, of course):

	30°C	60°C
pH 2.0		
pH 11.3		

Column Stability at pH 2.0



ACQUITY UPLC[™] BEH C18, 2.1 x 50 mm 30% ACN (v/v) with 0.1% TFA (pH 2.0) at 0.9 mL/min Analytes: protriptyline, amitriptyline, butyrophenone

Column Stability at pH=11.3

©2005 Waters Corporation

Waters



Acetonitrile-Pyrrolidine Buffer (pH=11.3) 45:55 (v/v) at 0.9 mL/min flow rate Analytes: butyrophenone, protriptyline and amitriptyline

©2005 Waters Corporation

How do we take advantage of smaller particles?

- Pressure tolerant particle/column
- High pressure fluidic modules
- Holistic design, low system volume, optimized flow paths
- Reduced cycle time autosampler with minimum carryover
- High speed detectors; optical and mass
- Controlled and coordinated system interaction
 - Software designed for system integration
 - Comprehensive diagnostic suite

ACQUITY UPLC[™] Technology Ultra performance by design

©2005 Waters Corporation

Detectors: Optical and/or Mass Spec Tunable UV or Photodiode Array Optimized flow cell for UPLC[™] High speed data sampling for

System Considerations: Small Footprint Redesigned tubing and fittings Consolidated waste management Integrated system diagnostics Connections Insight[™] remote diagnostics

Sample Organizer: (option) Expands capacity (22/15/8) Shuttles plate feed Heated/chilled

Sample Manager: Low dispersion XYZZ' Format Fast cycle times Low carryover Plates and/or vials Optional Sample Organizer Column Manager: Innovative pivot design Positions column to detector E-Cord[™] connection

Binary Solvent Manager:

High pressure blending Binary gradients Four solvent choices On-line degassing Low dispersion design UPLC pressure capabilities

©2005 Waters Corporation

ACQUITY UPLCTM Binary Solvent Manager

- Binary gradient
 - Independent style piston drives
 - Serial flow for each solvent, parallel gradient formation
 - High Pressure Gradient
 - Solvent select, 4 solvent choices
 - On board solvent degasser
 - 50µl mixer, ~120µL total system volume
 - Ability to deliver non-linear gradients
 - UPLC pressure capabilities
 - 0-1mL/min to 15K psi,
 - 1 -2mL/min to 9K psi
 - Dedicated transducer for each head.
 Enhances performance, provides diagnostic information



BSM Components



©2005 Waters Corporation

From Serial to Parallel Flow

- Two independently driven primary and accumulator piston drives
 - Separate drive motors
 - Two transducers
 - Two inlet check valves
- Four Solvents Capability
 - A1 or A2 and B1 or B2
 - Two solvents per gradient
- Six Channel Degasser
 - 4 chromatographic eluents
 - 2 solvents for Sample Manager
 - Strong & weak wash



Pressure Envelope

©2005 Waters Corporation

- Fluidic components optimized to support the pressure demands of UPLC[™] separations
 - Seals, head, bolts, check valves, fittings, tubing, mixer, pressure transducers
- UPLC[™] compatible pressure with flow rate range



©2005 Waters Corporation

- Glass Bead Mixing Technology
 - Tortured Mixing



Tee/Filter/Mixer

- Durable and will not disintegrate or affect delicate compounds
- Highly resistant to attack by most liquids and/or vapors.
- Glass Bead Mixer
 - Standard 4 x 5, approximate volume 50 μ L (P/N 700002631)
 - Optional Mixer 6 x 6, approximate volume 100 μ L (P/N 205000252)
- Filter
 - Similar to Alliance, 2 micron filter (P/N 700002674)



©2005 Waters Corporation

New Tubing & Fittings

- HPLC Fittings
 - Traditional drawn tubing



- Traditional Waters Fittings



- UPLC[™] Fittings
 - Center-less ground tubing



- Gold Compression Screw
- 2 Piece Swagelock Ferrules



ACQUITY 9/1000 SS316 Tubing

©2005 Waters Corporation




©2005 Waters Corporation

Standard 5/1000 inch 316 SS Tuning

@75 Times

@500 Times



New ACQUITY 2/1000 316 SS Tubing

©2005 Waters Corporation

@75 Times

@1000 Times



ACQUITY UPLCTM Binary Solvent Manager gradient steps

©2005 Waters Corporation

Naters

N = 6



ACQUITY UPLC[™] Binary Solvent Manager – Auto Blend

©2005 Waters Corporation



Injections of Caffeine with Dial-a-Mix at 1.0mL/min 4.0% to 5.0% ACN in 0.1% Increments

©2005 Waters Corporation

Original Assay Method from the USP for: Benzocaine, Butamben, and Tetracaine Hydrochloride Topical Solution

- *Mobile phase* Prepare a mixture of methanol, water, and 0.25 *M* sodium 1-heptanesulfonate (500:500:20). Make adjustments if necessary (see *System Suitability* under *Chromatography* (621).
- *Chromatographic system* (see *Chromatography* [621])—The liquid chromatograph is equipped with a 313-nm detector, and a 3.9-mm × 30-cm column that contains packing L1. The flow rate is about 2 mL per minute.
- Standard preparation— Transfer about 140 mg of <u>USP Benzocaine RS</u>, accurately weighed, to a 100-mL volumetric flask with the aid of 25 mL of methanol, and swirl. Transfer about 140 J mg of <u>USP Butamben RS</u>, accurately weighed, to the same volumetric flask with the aid of 25 mL of water, J being the ratio of the labeled amount, in percent, of butamben to the labeled amount, in percent, of benzocaine in the Topical Solution. Transfer about 140 J¢ mg of <u>USP Tetracaine Hydrochloride RS</u>, accurately weighed, to the same volumetric flask with the aid of 25 mL of water, J¢ being the ratio of the labeled amount, in percent, of benzocaine in the Topical Solution. Transfer about 140 J¢ mg of <u>USP Tetracaine Hydrochloride RS</u>, accurately weighed, to the same volumetric flask with the aid of 25 mL of water, J¢ being the ratio of the labeled amount, in percent, of tetracaine hydrochloride to the labeled amount, in percent, of benzocaine in the Topical Solution. Sonicate for about 1 minute, dilute with *Diluent* to volume, and mix.
- Procedure— Separately inject equal volumes (about 10 µL) of the Standard preparation and the Assay preparation into the chromatograph, record the chromatograms, and measure the peak areas for the major peaks.
- Chromatograph the Standard preparation, and record the peak responses as directed for Procedure: the relative
 retention times are about 0.3 for benzocaine, 0.8 for butamben, and 1.0 for tetracaine; the resolution, *R*, between
 the benzocaine peak and the butamben peak and between the butamben peak and the tetracaine peak is
 not less than 2; and the relative standard deviation for replicate injections is not more than 2.0% for each of
 the three analyte peaks.

ACQUITY UPLC[™] BSM Performance- RT Reproducibility

©2005 Waters Corpora



Component Summary Area RT Report

Reported by User: Andrew Aubin (Aubin)

Project Name: BBT

	SampleName	Inj	Vial	Benzocaine	Butamben	Tetracaine		
1	Ben Tetra Butarr	10	2:F,1	0.234	0.366	0.487		
2	Ben Tetra Butarr	11	2:F,1	0.236	0.371	0.493		
3	Ben Tetra Butarr	12	2:F,1	0.235	0.370	0.492		
4	Ben Tetra Butarr	13	2:F,1	0.235	0.369	0.492		
5	Ben Tetra Butarr	14	2:F,1	0.235	0.369	0.492		
6	Ben Tetra Butarr	15	2:F,1	0.235	0.370	0.492		
7	Ben Tetra Butarr	16	2:F,1	0.234	0.368	0.491		
8	Ben Tetra Butarr	17	2:F,1	0.234	0.366	0.487		
9	Ben Tetra Butarr	18	2:F,1	0.236	0.371	0.494		
10	Ben Tetra Butarr	19	2:F,1	0.235	0.369	0.492		
11	Ben Tetra Butarr	20	2:F,1	0.235	0.369	0.492		
Mean				0.235	0.369	0.491		
Std. Dev.				0.001	0.002	0.002		
% RSD				0.250	0.426	0.426		

Component Summary For Retention Time

Standard Deviations of 35, 94, 125 milliseconds

Alliance® HPLC Analytical Performance with isocratic conditions

©2005 Waters Corporation

Separation of Benzocaine, Butamben, and Tetracaine



ACQUITY UPLC[™]

Analytical Performance with isocratic conditions

©2005 Waters Corporation





Speed, Sensitivity, Resolution

©2005 Waters Corporation



©2005 Waters Corporation

Changes to Original Assay Method from the USP for: Benzocaine, Butamben, and Tetracaine Hydrochloride Topical Solution

- Mobile phase— Prepare a mixture of acetonitrile, water containing 10 mMol ammonium bicarbonate at a pH of 10.0 (1:1) for UPLC, 4:6 for HPLC. Make adjustments if necessary (see System Suitability under Chromatography (621).
- Chromatographic system (see Chromatography [621])—The liquid chromatograph is equipped with a 301-nm detector, and a 4.6-mm × 25-cm column that contains packing L1 (XTerra RP18, 5 um) for HPLC and a 2.1 X 50 mm 1.7 um for UPLC. The flow rate is about 1 mL per minute for UPLC and 2 mL per minute for HPLC.
- Standard preparation— Transfer about 140 mg of <u>USP Benzocaine RS</u>, accurately weighed, to a 100-mL volumetric flask with the aid of 25 mL of methanol, and swirl. Transfer about 140 J mg of <u>USP Butamben RS</u>, accurately weighed, to the same volumetric flask with the aid of 25 mL of water, J being the ratio of the labeled amount, in percent, of butamben to the labeled amount, in percent, of benzocaine in the Topical Solution. Transfer about 140 J¢ mg of <u>USP Tetracaine Hydrochloride RS</u>, accurately weighed, to the same volumetric flask with the aid of 25 mL of water, J¢ being the ratio of the labeled amount, in percent, of tetracaine hydrochloride to the labeled amount, in percent, of benzocaine in the Topical Solution. Transfer about 140 J¢ mg of USP Tetracaine Hydrochloride RS, accurately weighed, to the same volumetric flask with the aid of 25 mL of water, J¢ being the ratio of the labeled amount, in percent, of tetracaine hydrochloride to the labeled amount, in percent, of benzocaine in the Topical Solution. Sonicate for about 1 minute, dilute with Diluent to volume, and mix. Further dilute this stock 5 in100 for UPLC use.
- Procedure— Separately inject equal volumes (about 2 µL) of the Standard preparation and the Assay preparation into the chromatograph, record the chromatograms, and measure the peak areas for the major peaks.
- Chromatograph the *Standard preparation*, and record the peak responses as directed for *Procedure:* the relative retention times are about 0.3 for benzocaine, 0.8 for butamben, and 1.0 for tetracaine; the resolution, *R*, between the benzocaine peak and the butamben peak and between the butamben peak and the tetracaine peak is not less than 2; and the relative standard deviation for replicate injections is not more than 2.0% for each of the three analyte peaks.

ACQUITY UPLC[™]

Analytical Performance with Gradient conditions



Column: Mobile Phase A: Mobile Phase B: Flow Rate: Gradient: Wash Solvents:

Wavelength: Data Rate: 2.1 x 50 mm @ 65° C 0.1% Formic Acid in Water 0.1% Formic Acid in ACN 1.3 mL/min 30% - 90% B over 0.5minutes, Curve 5 Strong – 50/50 Water/ACN (200µL); Weak – 95/5 Water/ACN 600µL) 240nm 40 pts/s, Filter Constant: 0.0

Average T _R %RSD	0.109%
Average Peak Area %RSD	0.368%
Average 4 Sigma Peak Width	0.369 sec
Average w _{1/2} Peak Width	0.210 sec
Minimum Resolution	2.15
Total Cycle Time	56 sec
Maximum Pressure	13,500 psi

ACQUITY UPLCTM Analytical Performance with Gradient conditions



©2005 Waters Corporation

ACQUITY UPLC[™] Non-Linear Gradient Capability

Gradient Profiles on ACUITY UPLC™



©2005 Waters Corporation

Method Development Challenge:

Gradient elution of 9 component mix

- To go as fast as possible
- To meet a minimum resolution requirement of 2.0
- Start with generic method 30s gradient with 15s hold at final conditions. (5% to 95% ACN)
- Run gradient using curves 3-9

Chromatograms of Curves 3-9

©2005 Waters Corporation

Waters

©2005 Waters Corporation

Curve	3	4	5	6	7	8	9
Resolution	1.32	1.76	2.17	3.01	3.59	3.61	3.36
Elution Time of Last Peak	0.34	0.40	0.46	0.57	0.64	0.66	0.68

Shortest Run Time with required resolution is using Curve 5

©2005 Waters Corporation

ACQUITY UPLC[™] Sample Manager

New Technologies - *all-new* injection process

- Pressure tolerant, low volume flowthrough isolation & injection valves
- Fast cycle times, low carry over
- 2 Plate capacity, multiple formats
 Well plates, microtube plates, vials
- Thermal control (4 40°C)
- Low dispersion column management
- Sample Organizer (optional)

©2005 Waters Corporation

"Pod" Style Injector

• Modular design (removable "Pod") permits easy access

ACQUITY UPLC[™] Loops available

©2005 Waters Corporation

2 μL 5 μL 10 μL **20 μL** 50 μL

ACQUITY UPLC[™] Sample Manager Performance - Precision

©2005 Waters Corporation

Component Summary Area RT Report

Reported by User: Andrew Aubin (Aubin)

Project Name: PITTCON_2005

Component Summary For Area

	SampleName	Date Acquired	lnj Vol (ul)	Inj	Vial	3,4-MDA	2,3-MDA	3,4-MDMA	2,3-MDMA
1	2,3+3,4 MDA MDMA	2/14/2005 10:28:57 AM	5.00	1	2:D,1	252746	169926	235323	164552
2	2,3+3,4 MDA MDMA	2/14/2005 10:33:42 AM	5.00	1	2:D,1	251676	169802	235335	163520
3	2,3+3,4 MDA MDMA	2/14/2005 10:36:38 AM	5.00	1	2:D,1	252557	169395	235366	163371
4	2,3+3,4 MDA MDMA	2/14/2005 10:39:35 AM	5.00	1	2:D,1	251382	169090	232696	162642
5	2,3+3,4 MDA MDMA	2/14/2005 10:42:31 AM	5.00	1	2:D,1	251354	168659	234721	163114
6	2,3+3,4 MDA MDMA	2/14/2005 10:45:27 AM	5.00	1	2:D,1	250666	169286	234564	163830
7	2,3+3,4 MDA MDMA	2/14/2005 10:48:24 AM	5.00	1	2:D,1	252694	168439	235468	163120
8	2,3+3,4 MDA MDMA	2/14/2005 10:51:20 AM	5.00	1	2:D,1	252331	169251	235365	163195
9	2,3+3,4 MDA MDMA	2/14/2005 10:54:14 AM	5.00	1	2:D,1	252360	168626	234133	164507
10	2,3+3,4 MDA MDMA	2/14/2005 10:57:11 AM	5.00	1	2:D,1	252659	168037	234904	162978
Mean						252043	169051	234787	163483
Std. Dev.						722	604	856	636
% RSD						0.29	0.36	0.36	0.39

ACQUITY UPLC[™] Sample Manager Performance - Linearity

©2005 Waters Corporation

0.1 mg/mL Acetaminophen with 20µL loop in Pressure-Assist Mode

ACQUITY UPLC[™] Sample Manager Performance -Carryover

©2005 Waters Corporation

Wash Conditions: 200µL Strong (6:3:1 / 0.1% phosphoric in Water:ACN:IPA) 600µL Weak (0.01M HCI)

$\begin{array}{l} ACQUITY \ UPLC^{\mathsf{TM}} \\ Sample \ Manager \ Performance \ -Carryover \end{array}$

©2005 Waters Corporation

Wash Conditions: 200µL Strong (6:3:1 / 0.1% phosphoric in Water:ACN:IPA) 600µL Weak (0.01M HCI)

©2005 Waters Corporation

Cycle Times?

ACQUITY UPLC™ Injection-to-Injection Cycle Times

Injection-to-injection cycle times are calculated as the time it takes for the injector to perform all of the tasks during the injection cycle from the point of initiation by the data acquisition software until the valve turns and the injection is made. To determine this value a series of injections are monitored by a detector that is continually collecting data, even between the run times. The difference in time between corresponding peaks of each injection gives the total cycle time. Subtracting the run time from the total cycle time will yield the injection-to-injection cycle time.

	8 8		2 p
Injection Interval	Total Cycle Time		Injector Cycle Time
0.65 8		(Se©)	
0.60-	64.9	[#] 36.0	28.9
0.55- 0.50- 2	65.1	36.0	29.1
0.45-3	65.0	36.0	29.0
0.40 ح 0.35	64.1	36.0	28.1
0.00 5	63.0	36.0	27.0
0.25 6	65.5	36.0	29.5
7182.11 2220-11	6 % 0. %	36 0	⁸⁸ 31∮0
0.10-	66.9	BF 36.0	30.9
-0.05 Average	65.2	36.0	29.2
0.00 1.00 2.00	300 400 500	600 Z00 800	900 1000 1100 1200

©2005 Waters Corporation

ACQUITY UPLCTM Column Manager

Column Manager

- Symbiotic with Sample Manager
 - Flow path distance optimized
 - 65°C upper limit
 - Up to 150 mm column
- Pivot positioning
 - Stacked mode
 - w/ optical detector
 - Swung out mode
 - w/ MS detector only
 - Column may be accessed from either side

©2005 Waters Corporation

Attaching The Column to the System

High Pressure Finger Tight Fittings Holds >15000 psi

Patent Pending

Column Outlet

©2005 Waters Corporation

- One piece finger tight fittings work well
- Cut 0.0025" Peek to length
- Keep as short as possible to minimize dispersion

©2005 Waters Corporation

- Column heating considerations
 - Consistent temperature for RT stability
 - Solvent viscosity is lower at higher temperatures
 - Lower backpressures at higher temperatures

ACQUITY UPLC[™] Column Manager with eCord[™]

©2005 Waters Corporation

$\mathbf{eCord^{\mathsf{TM}}} \ \mathbf{Technology}$

- Column Information

 Certificate of Analysis
- Column Usage History
 - Total number of samples
 - Injections
 - Pressure/Temperature history
 - And much more...

eCord[™]

©2005 Waters Corporation

ACQUITY UPLC[™] eCord[™]

Ultra Performance LC Columns

TM

eCord™ Column QC

©2005 Waters Corporation

- Faster eluting peaks require higher data rates
 - High light throughput/transmission
 - Faster digital filter time constants
- UPLC[™] separations require low dispersion flow cells
 - Maintain chromatographic fidelity
- Tunable UV, PDA and ELSD Options
- MS Options

ACQUITY UPLC[™]

Detectors

©2005 Waters Corporation

ACQUITY UPLCTM Flow Cell Design

- Light Guided UPLC[™] flow cells
 - 10 mm pathlength, 500 nL volume
 - Flow cell channel is the inside of a low-index Teflon AF tube
 - Total internal reflection at walls, like optical fiber cladding

ACQUITY UPLC[™] Low dispersion UPLC optics

©2005 Waters Corporation

Low dispersion UPLCTM optics Expanded Region

©2005 Waters Corporation

Data Acquisition Rates Impact on UV chromatography data



Effect of Filtering Constant on Chromatography

©2005 Waters Corporation



5 component Analgesic Mix on 2.1x50mm ACQUITY UPLC[™] Column, 30s gradient

UV Spectra on ACQUITY PDA



Peak Purity and Library Match Capabilities



©2005 Waters Corporation

ACQUITY UPLC[™] Optional Sample Organizer

Supports high throughput/high capacity

- Expands capacity-mix or match
 - 22 thin (microtitre) plates
 - 8,448 samples max
 - 15 medium height plates
 - 8 full height (vial) plates
 - 384 x 2mL vials
- Temperature controlled
- Programmable Access to any sample location
- Supports Walk-up Operation



ACQUITY UPLC[™] Optional FlexCart

©2005 Waters Corporation



Supports inlet mobility

©2005 Waters Corporation

How do we take advantage of smaller particles?

- Pressure tolerant particle/column
- High pressure fluidic modules
- Holistic design, low system volume, optimized flow paths
- Reduced cycle time autosampler with minimum carryover
- High speed detectors; optical and mass
- Controlled and coordinated system interaction
 - Software designed for system integration
 - Comprehensive diagnostic suite

System Information



Sample Syringe Diagnostic



Needle Seal - Static vs Dynamic TEST



©2005 Waters Corporation

ACQUITY UPLCTM Designed with Intent

UPLC[™] enabled, laboratory rugged

- Performance and Reliability
 - UPLC pressure capabilities
 - Low dispersion fluidics
 - Rapid, clean injections
 - Fast detection
 - Rugged chemistry
 - Analytical performance
 - Extensive onboard diagnostics
 - Connections Insight[™]

- Convenience
 - AQT Qualification
 - Small footprint
 - Integrated waste management
 - eCord[™] Technology
 - Intuitive Control Console
 - Multiple plate formats
 - High capacity sample organizer
 - FlexCart

UPLC[™] = a higher quality of (more) information faster



©2005 Waters Corporation

QUESTIONS?

TM COULT Ultra Performance LC







©2005 Waters Corporation

Acknowledgements

- John Morawski
- Tanya Jenkins
- Eric Grumbach
- Robert Plumb
- Jeffrey Mazzeo