

# 13 Introduction to Inlets

## **Inlet types**

### **Using hydrogen**

Procedure: Pressure units: Select psi, kPa, bar

### **The inlet and column control tables**

#### **The column control tables**

The column control table—defined capillary columns

The column control table—packed or undefined capillary columns

#### **What is gas saver?**

Procedure: Using gas saver

#### **Pre Run and Prep Run**

The [Prep Run] key

Procedure: Auto Prep Run

#### **Septum purge**

# Introduction to Inlets

---

## Inlet types

The 6890 GC has five types of inlets available. All are offered with electronic pneumatics control (EPC) and two are offered without.

**Table 28 Inlet Types**

Inlet type	Gas control
Split/splitless	EPC and nonEPC
Purged packed	EPC and nonEPC
Cool on-column	EPC only
Programmed temperature vaporization	EPC only
Volatiles interface	EPC only

---

## Using hydrogen

---

### **WARNING**

When using hydrogen (H<sub>2</sub>), as the carrier gas, be aware that hydrogen (H<sub>2</sub>) gas can flow into the oven and create an explosion hazard. Therefore, be sure that the supply is off until all connections are made, and ensure that the inlet and detector column fittings are either connected to a column or capped at all times when hydrogen (H<sub>2</sub>) gas is supplied to the instrument.

---

---

### **WARNING**

Hydrogen (H<sub>2</sub>) is flammable. Leaks, when confined in an enclosed space, may create a fire or explosion hazard. In any application using hydrogen (H<sub>2</sub>), leak test all connections, lines, and valves before operating the instrument. Always turn off the hydrogen (H<sub>2</sub>) supply at its source before working on the instrument.

---

**Table 29 An Overview of Inlets**

<b>Inlet</b>	<b>Column</b>	<b>Mode</b>	<b>Sample concentration</b>	<b>Comments</b>	<b>Sample to column</b>
Split/splitless	Capillary	Split Pulsed split	High High	May be useful with large injections	Very little Very little
		Splitless Pulsed splitless	Low Low	Useful with large injections	All All
Cool on-column	Capillary	n/a	Low or labile	Minimal discrimination and decomposition	All
Purged packed	Packed Large capillary	n/a	Any	OK if resolution not critical	All
		n/a	Any		All
Programmed temperature vaporization	Capillary	Split	High	Multiple injections concentrate analytes and vent solvent	Very little
		Pulsed split	High		Very little
		Splitless	Low		All
		Pulsed splitless	Low		All
		Solvent vent	Low		Most
Volatiles interface	Capillary	Direct	Low	Lowest dead volume	All
		Split	High	Max flow = 100 mL/min	Very little
		Splitless	Low		All

**Table 30 Column Size and Carrier Gas Flow Rate**

Column type	Column size	Carrier gas flow rate	
		Hydrogen	Helium
Packed	1/8-inch		30
	1/4-inch		60
Capillary	50 $\mu\text{m}$ id	0.5	0.4
	100 $\mu\text{m}$ id	1.0	0.8
	200 $\mu\text{m}$ id	2.0	1.6
	250 $\mu\text{m}$ id	2.5	2.0
	320 $\mu\text{m}$ id	3.2	2.6
	530 $\mu\text{m}$ id	5.3	4.2

These flow rates, in mL/min at normal temperature and pressure (25°C and 1 atm) are recommended for all column temperatures.

For capillary columns, flow rates are proportional to column diameter and are 20% lower for helium than for hydrogen.

### Procedure: Pressure units: Select psi, kPa, bar

You can display pressure in psi, bar, or kPa. To check the units you are using, pressing the [Info] key while the cursor is on the Pressure line of a control table. To change the display units:

1. Press [Options].
2. Scroll to **Keyboard & Display** and press [Enter].

```

      OPTIONS
  Calibration
  Communication
  Keyboard & Display <
  Diagnostics
  
```

3. Scroll to **Pressure units:** and press [Mode/Type].

```

Keyboard lock      Off
Key click         On
Warning beep      On
Method mod beep   Off
  KEYBOARD OPTIONS
  Pressure units:  psi <
  Radix type:     .
  
```

```

      PRESSURE UNITS
  *psi
  bar
  kPa
  
```

4. Choose a new pressure unit and press [Enter].

**Table 31 Pressure Unit Conversions**

To convert	to	Multiply by
psi	bar	0.0689476
	kPa	6.89476
bar	psi	14.5038
	kPa	100
kPa	psi	0.145038
	bar	0.01

---

## The inlet and column control tables

The tables for the inlet and column are interrelated. If you set a pressure at the column control table, that same pressure setting is active on the inlet control table, and vice versa. Although pneumatics can be controlled from either the column or the inlet, the column should be considered first.

COLUMN 1 (He)	
Dim	30.0 m 320 u
Pressure	10.0 10.0
Flow	0.7
Velocity	19
Mode:	Constant flow

FRONT INLET (S/SL)	
Mode:	Splitless
Temp	250 250 <
Pressure	10.0 10.0
Purge time	0.75
Purge flow	15
Total flow	??
Gas saver	Off

Note that the pressure readings—both setpoint and actual—are identical on the column and inlet control tables.

## The column control tables

The control tables change depending on your column configuration. The next few pages describe the column control tables for the two types of columns, capillary and packed.

### The column control table—defined capillary columns

If your column is defined, your control table will be similar to [Figure 38](#).

*The title* This heading identifies the column—Column 1 or Column 2— and the type of carrier gas configured to the inlet (in parentheses).

*Dim* This line shows the column dimensions you have specified. Column length is in meters (m) and column inside diameter is in microns ( $\mu$ ).

Pressure, flow, and velocity are related. If the column is defined, enter any one of them and the GC computes and displays the other two.

*Pressure* The setpoint appears at the far right. The number at the left shows the actual pressure value. When you enter a pressure value, the values for flow and average linear velocity are calculated and displayed.


*Flow* If you enter a flow (in mL/min) here, pressure and velocity are calculated and adjusted.

*Velocity* If you enter average linear velocity (in cm/sec), pressure and flow are calculated.

*Mode:* There are four column modes: constant flow, constant pressure, ramped flow, and ramped pressure. To change the mode, scroll to *Mode:* and press [Mode/Type].

[“Flow and Pressure Control”](#) explains how to set pressure and flow programs.

Press [Col 1] or [Col 2]



**Mode:** Your control table also has one of these, depending on Mode:

Mode: Const flow	<
------------------	---

Mode: Const pressure	<
----------------------	---

Mode:Ramped flow	<
Init flow	4.0
Init time	2.0
Rate 1	0.5
Final flow 1	8.0
Final time 1	2.0
Rate 2 (Off)	0.00

Mode:Ramped pressure	<
Init pressure	10.0
Init time	1.0
Rate 1	1.0
Final pressure	125.0
Final time 1	5.0
Rate 2 (Off)	0.00

**Figure 38 Column display — defined capillary columns**

## The column control table—packed or undefined capillary columns

If you have not defined your column or if your inlet selection is *Unspecified*, your column control table will be similar to [Figure 39](#).

*The title* This heading identifies the column—Column 1 or Column 2— and the type of carrier gas configured to the inlet (in parentheses).

*Dimensions unknown* This line tells you that you have not defined your column.


*Pressure* The *split/splitless* inlet and the *cool on-column* inlet are pressure controlled. Because the column is unknown, flow and average linear velocity cannot be computed.

The *purged packed* inlet is flow controlled. The actual pressure is displayed, but is not controllable by the user.


*Mode:* You have a choice of three modes if using a *split/splitless* or *cool on-column* inlet—constant pressure, constant flow, and ramped flow. The *packed* inlet gives you only the two flow modes—constant and ramped.

[“Flow and Pressure Control”](#) explains how to set pressure and flow programs.

**Split/splitless or cool on-column inlets**



**Purged packed inlet**




**Figure 39 Column display — Packed or undefined capillary columns**

## What is gas saver?

Gas saver reduces carrier gas flow from the split vent after the sample is on the column. Column head pressure and flow rate are maintained, while purge and split vent flows decrease. Flows—except column flow—remain at the reduced level until you press [Prep Run].

You can use gas saver in all modes of operation of the Split/Splitless and PTV inlets and in the split and splitless modes of the Volatiles Interface.



**Figure 40 Gas saver operation**

The pulsed modes of the split/splitless and PTV inlets are similar except for the pressure pulse starting at [Prep Run] and ending at Pulse time. The solvent vent mode of the PTV is more complex; see [“Using the Solvent Vent Mode”](#) for details.

### Procedure: Using gas saver

Press [Front Inlet] or [Back Inlet].

Mode:		Split
Temp	24	Off
Pressure	0.0	Off
Split ratio		10
Split flow		0.0
Tot flow	0.0	Off
FRONT INLET (S/SL)		
Gas saver		On
Saver flow		20.0
Saver time		2.00

1. Turn on gas saver.

2. Set a flow. Must be at least 15 mL/min greater than the column flow.

3. If in split mode, set after injection time.  
In all other modes, set after purge time.

## Pre Run and Prep Run

With some inlets and operating modes, certain instrument setpoints are different between runs than during an analysis. To restore the setpoints for injection, you must place the GC into the Pre Run state.

You must use the Pre Run state when:

- Using gas saver with any inlet.
- Using splitless mode with any inlet.
- Using a pressure pulse mode with any inlet.
- Using the solvent vent mode of the PTV inlet.
- Using the direct or splitless mode of the Volatiles Interface.

There are two ways to begin Pre Run—manually push the [Prep Run] key before each run or configure the GC to enter the Pre Run state automatically. The two methods are discussed below.

During the Pre Run state:

- The Pre Run light blinks and Not Ready is on.
- Setpoints change to the correct values for injection.
- Inlet, detector, and oven equilibration times begin.

When all equilibration times expire, the Pre Run light is on steadily. When all criteria for a run are met, the Not Ready light turns off. The GC is now ready for sample injection.

### The [Prep Run] key

Press the [Prep Run] key before you inject a sample manually. The GC enters the Pre Run state. When the Pre Run light is steady and the Not Ready light goes off, begin the analysis.

### Procedure: Auto Prep Run

With most automatic injection systems, you do not need to use the [Prep Run] key. If your sampler or automation controller (for example, an integrator or workstation) does not support the [Prep Run] function, you must set the GC to Auto Prep Run. To do this:

1. Press the [Config] key to view a list of configurable parameters.

2. Scroll to the Instrument parameter and press [Enter].
3. Scroll to Auto prep run and press [On].

```

CONFIG INSTRUMENT
Serial#US00100001
Auto prep run      On <
F inlet type      None
B inlet type      PP
  
```

---


## Septum purge

The septum purge line is near the septum where the sample is injected. A small amount of carrier gas exits through this line to sweep out any bleed.

Each inlet has a different septum purge flow. The GC automatically sets the purge flow for EPC inlets, but you can measure it from the septum purge vent at the flow manifold if you like.

**Table 32 Septum Purge Flows**

Inlet	Carrier	Septum purge (mL/min)
Split/splitless, all modes	He, N <sub>2</sub> , Ar/5%Me	3
	H <sub>2</sub>	6
Purged packed	All	1 to 3
Cool on-column	He, N <sub>2</sub> , Ar/5%Me	15
	H <sub>2</sub>	30
PTV	He, N <sub>2</sub> , Ar/5% Me	3
	H <sub>2</sub>	6
Volatiles interface	He, N <sub>2</sub> , Ar/5%Me	3
	H <sub>2</sub>	6



**Figure 41** Septum purge vents