

Appendix C

The MassLynx Tune Window

Note: For the highest mass accuracy, tune and calibrate the instrument using a suitable reference compound before sample data are acquired (see Chapter 2, Chapter 3 and Chapter 6.)

C.1 Opening the MassLynx Tune Window

Click the MassLynx window Instrument shortcut bar MS Tune icon to open the Tune window. Figure C-2 summarizes the Tune window layout.

C.2 Selecting the Ionization Mode

Select the required ionization mode from the Ion Mode menu (Figure C-1). The current mode has a check mark next to it.




Figure C-1 Ion Mode Menu

The screenshot shows the MassLynx Tune Window interface. At the top is a menu bar with options: File, Ion Mode, Calibration, Gas, Ramps, Options, Help. Below the menu bar is a toolbar with icons for saving, opening, printing, and displaying information. A 'Toolbox' callout explains that some routine operations can be performed with a single click. A 'Current Gain' gauge is located in the top right corner, showing a value of 1.00e-4.

The main window is divided into several sections:


- Parameters Section:** Includes fields for 'CapRay (V)', 'Cone (V)', 'Extractor (V)', 'RF Lens (V)', 'Temperatures' (Source Temp, Desolvation Temp), 'Gas Flow' (Desolvation Flow), and 'Syringe' (Pump Flow, Syringe Status).
- Function Table:** A table with columns for 'Function', 'Set', 'Mass', 'Scan', and 'Gain'. It lists functions like 'MS Scan' and 'Trap Scan'.
- Peak Display:** A mass spectrum plot showing intensity versus mass-to-charge ratio (m/z). The x-axis ranges from 453.0 to 457.0. A 'Peak Display' callout indicates that up to four masses can be displayed and zoomed.
- Callouts:** Several callouts provide detailed instructions:
 - 'Click to display the tune parameters for that region' points to the parameter fields.
 - 'Check up to four boxes to display the peaks' points to the function table.
 - 'Click on the arrow to select the scan function' points to a dropdown arrow in the function table.
 - 'Slider Bar' points to a control for the 'Desolvation Flow'.
 - 'Readback Window' points to a window showing the current value of the 'Desolvation Flow'.
 - 'Display tune peak information' points to a toolbar icon.
 - 'Display vacuum information' points to another toolbar icon.
 - 'Toggle on/off collision energy ramp' points to a toolbar icon.
 - 'Print current window in portrait format' points to a toolbar icon.
 - 'Open an existing tune file' points to a toolbar icon.
 - 'Enabled for MS-MS functions' points to a checkbox in the function table.

At the bottom of the window, there is a status bar with labels for 'Results', 'Vacuum OK', and 'Currents'.


Figure C-2 Tune Window

C.3 Controlling Gas Flows

C.3.1 Controlling the Nebulizer, Desolvation, and Cone Gas Flows

To toggle the nebulizer, desolvation, and cone gasses on and off, click  or select Gas > Gas.

C.3.2 Controlling the Collision Gas Flow

To toggle the collision gas on and off, click  or select Gas > Collision Gas.

C.4 Controlling the Syringe Pump

To toggle the syringe pump on and off, click .

C.5 Selecting the Syringe Type

1. Select Options > Syringe Type to open the Syringe Selection dialog box (Figure C-3).





Figure C-3 Syringe Selection Dialog Box

2. Choose the syringe type from the drop-down list.
If the required syringe type is not available in the drop-down list, click Edit to open the syringe list, and then add the required details.
3. Click OK.

C.6 Selecting the Scan Time and Inter Scan Delay

Enter the Scan Time and Inter Scan Delay in the Scope Setup dialog box (Figure C-4).

This is opened by clicking  or selecting Options > Scope Parameters.




Figure C-4 Scope Setup Dialog Box

The Scan Time(s) and Inter Scan Delay (s) parameter values control the speed with which the tune peak display is updated.

Tuning is more responsive when these parameters are set to low values.

C.7 Setting the Ramp Controls

C.7.1 Creating a Cone Voltage Ramp

1. Select Ramps > Cone Ramp Gradient to open the Cone Ramp dialog box (Figure C-5).

Two values of cone voltage (Cone Start Volts and Cone End Volts) are defined at two particular masses (Start Mass and End Mass). These values define a gradient for the cone voltage, which is then extrapolated to cover the full mass range.

2. Make any changes required, and click OK to exit.





Figure C-5 Cone Ramp Dialog Box

C.7.2 Controlling the Cone Voltage Ramp

To toggle the cone voltage ramp on and off, click  or select Ramps > Use Cone Ramp.

C.7.3 Creating a Collision Energy Ramp

1. Select Ramps > Collision Energy Ramp Gradient to open the Collision Ramp dialog box (Figure C-6).





Figure C-6 Collision Ramp Dialog Box


Two values of collision energy (CE Start and CE End) are defined at two particular masses (Start Mass and End Mass). These values define a gradient for the collision energy voltage, which is then extrapolated to cover the full mass range.

2. Make any changes required, and click OK to exit.

C.7.4 Controlling the Collision Energy Ramp

To toggle the collision energy ramp on and off, click  or select Ramps > Use Collision Energy Ramp.

C.8 Resetting the Zero Level

The zero level (or baseline) can be repositioned by clicking  or by selecting Options > Reinitialize.

This command causes the instrument control system to measure the position of the noise signal so that any baseline offset caused by the electronics or instrumentation can be compensated for.

You should reset the zero level whenever the multiplier voltage is changed.

C.9 Controlling the Display of Readback Windows

The Readbacks dialog box (Figure C-7) is opened by selecting Options > Readbacks.




Figure C-7 Readbacks Dialog Box

There are three options for displaying system readbacks in the Tune window:

- Readbacks displayed continuously (Always On).
- Readbacks hidden (Always Off).
- Readbacks displayed only when differing from their defined values by more than 10% (On out of range).

A number of the readbacks are for indication purposes only and are not true (calibrated) records of the actual voltages on the instrument. The acceptable variation between the set value and the readback value varies depending on the particular tune parameter. If you are concerned about any reading, contact Waters for advice (see Section 8.6).

C.10 Changing Tune Parameter Settings

You can modify most parameters in the following ways:

- Use the mouse to drag the slider bar.
- Click the slider bar and use the left and right arrow keys to change the value by one increment. The edit window is updated as the slider bar is activated.
- Type a new value into the edit window, then press the Enter key.

Other parameters have only an edit window and are changed by direct typing.


The speed with which the system responds to changes depends on the speed with which the peak display is refreshed. For the fastest response, set the Scope Setup dialog box's Scan Time (s) and Inter Scan Delay (s) values as short as possible (see Figure C-4).

C.11 Saving Instrument Tune Parameters

Instrument tuning parameters can be saved in an instrument parameter file (.ipr), which can be recalled later.

An instrument parameter file contains all the parameters for all supported ionization modes, not just the ionization mode currently selected. Instrument parameter files also contain settings for the analyzer, inlet set points, and peak display.

C.11.1 Creating a New Instrument Parameter File

To create a new instrument parameter file, click  or select File > New.

C.11.2 Saving Tune Parameters in an Instrument Parameter File

To save the current tune parameters with the existing instrument parameter file name, click


 or select File > Save.

To save the current tune parameters with a new file name:


1. Select File > Save As to open the Save As dialog box.
2. Enter a new file name or select an existing file from the displayed list.
3. Click Save.
4. If the selected file already exists on disk, a warning is displayed. Click Yes to overwrite the existing information, or No to enter a different file name.

C.11.3 Opening an Existing Instrument Parameter File

To open an existing instrument parameter file, i.e., restore a saved set of tune parameters:

1. Click  or select File > Open to open the Open dialog box.
2. Select the required instrument parameter file, either by typing its name or selecting one from the list.
3. Click Open.

C.12 Printing Tune Information

To print a report containing a copy of the on-screen tune peak information and a record of each parameter setting, click  or select File > Print.

This report is not configurable by the user.

C.13 Using the EasyTune Source Page

The Tune window EasyTune Source page replaces the Source and Analyser pages. It displays the tune parameters commonly used during instrument tuning, while the more rarely used parameters are hidden.



By default, the MassLynx Tune window opens showing the Source and Analyser pages. To open the Tune Window EasyTune Source page (Figure C-8), select Options > EasyTune.

Using the EasyTune Source page, the instrument is tuned in a similar manner to that described in Chapters 2 and 3.

To close the EasyTune Source page and return to the Source and Analyser pages, select Options > Advanced.

C.14 Changing Between the Peak and Vacuum Displays

Peak information, or vacuum gauge readbacks, can be displayed in the top-right corner of the Tune window.

- To display the vacuum information, select  or select Options > Vacuum Monitor.
- To display the peak information, select  or Options > Peak Editor.

C.15 Setting-Up Peaks for Tuning

C.15.1 Selecting Peaks for Tuning

1. With the Tune window set to display peak information (see Section C.14), select the peaks to be displayed by selecting the appropriate numbered box(es) in the top-right of the Tune window.
2. For each active peak, enter the required values for the Mass, Span, and Gain.

Note: For MS/MS functions, Set is also enabled, allowing the mass of the parent, daughter, neutral loss, or neutral gain ion (as appropriate) to be entered.




Figure C-8 Tune Window EasyTune Source Page


C.15.2 Selecting the Operating Mode for a Peak

To select the operating mode (scan function), select the required mode for the peak from the Function drop-down list.

C.15.3 Selecting the Tune Mass for a Peak


Either:

1. In the appropriate peak window, on the right side of the Tune window, click and drag the mouse within the bounds of the axis to draw a "rubber band" around the region of interest.


2. Release the mouse button. This range is redisplayed to fill the window. The value displayed in the Mass box is the mass at the center of the window. This operation can be repeated as often as required.
3. Clicking  once displays the previous magnification range and mass; clicking it again returns to the default settings.


Or:

1. Enter a value in the Mass box for the required peak.
2. Click Return.

This becomes the default, so if the range is then altered using the mouse and you click  twice, Mass returns to this value.

Or:


1. Position the cursor at the top of the peak window, just below the line displaying the gain value (e.g., 2×).
2. When  appears, click and drag the mouse until the required mass is displayed in the Mass box and at the top of the window.

This becomes the default, so if the range is altered using the mouse and you click  twice, Mass returns to this value.

C.15.4 Selecting the Span of a Displayed Peak

Either:


1. Click and drag the mouse, horizontally, from one end to the other of the region of interest. As the mouse is dragged, a “rubber band” indicates the selected range. Do not go beyond the bounds of the axis.
2. Release the mouse button to redisplay the selected range so that it fills the current window. This operation can be repeated as often as required.

Clicking  once displays the previous magnification range, clicking it again returns to the default settings.

Or:

1. Enter a value in the Span box for the required peak.

2. Click Return.

This becomes the default, so if the range is altered with the mouse and you click  twice, Span returns to this value.

C.15.5 Changing the Gain of a Displayed Peak

Either:

1. Double-click the line displaying the gain value (e.g., 2×) above the peak, to double the gain applied to that peak.
2. Double-click below the peak display to halve the gain.

Or:

1. Click and drag the left mouse button, vertically, from one end to the other of the region of interest.

As the mouse is dragged, a marquee indicates the selected range.

Do not go beyond the bounds of the axis.

2. Release the mouse button to redisplay the selected range so that it fills the current window.

Or:

1. Enter a value in the Gain box for the required peak.
2. Click Return.

C.16 Customizing the Peak Display

C.16.1 Opening the Peak Display Menu

Use the Peak Display menu (Figure C-9) to customize the Peak Display. Open it by right-clicking in the appropriate peak display window in the Tune window.

Note: *The display window for each peak can be individually customized, e.g., the peak color for peak 1 can be red, for peak 2 green, etc.*




Figure C-9 Peak Display Menu

C.16.2 Customizing the Colors and Numbers of Displayed Traces

To change the color of the background and traces, and to change the number of traces displayed, select **Customise > Plot Appearance** from the Peak Display pop-up menu to open the Customise Plot Appearance dialog box (Figure C-10).




Figure C-10 Customise Plot Appearance Dialog Box

To change the colors on the display, click the color box adjacent to Newest Trace, Background, or Trace Fill as required, and select a new color from the Color dialog box.

To change the number of displayed traces, enter the required value in the Visible traces box, within the range 2 to 20.

If more than one trace is displayed, the older traces can be displayed in different color shades to the newer ones. Drag the Colour Interpolation slider toward the full position. The colors of the older traces appear in the Trace colour sample (new->old) field.

C.16.3 Customizing the Peak Trace Line Appearance

Each trace may be displayed as:

- Outline only – select Trace > Outline from the Peak Display pop-up menu.
- With the area below the line filled – select Trace > Fill from the Peak Display pop-up menu.
- Maximum and minimum points only – select Trace > Min/Max from the Peak Display pop-up menu.

The selected option displays a check mark in the Peak Display pop-up menu.

C.16.4 Customizing the Peak Intensity Display

To display the peak intensities as absolute values (counts/second), select Intensity > Absolute Intensity from the Peak Display pop-up menu.

To display the peak intensities as percentage values relative to the intensity of the highest peak, select Intensity > Relative Intensity from the Peak Display pop-up menu.

The Peak Display pop-up menu Intensity > Normalise Data option can be selected in conjunction with either of the above options. It controls the way in which the peak display is scaled. When enabled, the display scales to the value of the intensity of the highest peak; when disabled, the display scales to the default value set in the MassLynx software. It is recommended that the Intensity > Normalise Data option is normally enabled.

The selected options display a check mark in the Peak Display pop-up menu.

C.16.5 Customizing the Peak Display Grid

The Peak Display vertical and horizontal grid lines may be independently displayed or hidden.

To display the horizontal grid lines, select Grid > Horizontal from the Peak Display pop-up menu.

To display the vertical grid lines, select Grid > Vertical from the Peak Display pop-up menu.

The selected options display a check mark in the Peak Display pop-up menu.

C.17 Selecting the Instrument Name

1. Select Options > Instrument Name to open the Instrument Name dialog box.
2. Enter the required name in the ID text box.
3. Click OK.

C.18 The Diagnostics Page



Caution: *The instrument may be damaged if the parameters in the Tune window Diagnostics page are modified by an unqualified user. The parameters on this page should only be modified by a Waters field service engineer.*

The Tune window Diagnostics page, opened by selecting Options > View Diagnostics Page, is intended for use by Waters field service engineers only.

C.19 Manually Controlling the T-Wave Optics

The T-Wave devices have been optimized for use under normal operating conditions and their control should normally be left to MassLynx. Nevertheless, manual control of the optics may be desirable under certain conditions, for example, at high gas cell pressures ($>7 \times 10^{-3}$ mbar).

Note: *Inappropriate settings will result in poor performance in terms of sensitivity and/or high levels of cross-talk.*

The source T-Wave is only applied when using the optional MUX-technology interface (see the *Waters Micromass Quattro Premier XE MUX-technology Interface Operator's Guide* for details). The collision cell T-Wave is applied for MS methods that require collision gas in the cell.

To select manual control of the T-Wave devices, select Options > View T-Wave Manual Controls in the Tune window. This opens the T-WAVE page (Figure C-11). Separate control of the source and/or collision cell T-Wave pulse voltage and velocity is enabled by selecting the appropriate Enable Manual Controls check box(es).




Figure C-11 Tune Window T-WAVE Page

Each T-Wave device has two adjustable parameters: the Pulse Velocity (m/s) and Pulse Voltage. The standard setting for the Pulse Velocity is 300 m/s. Setting a velocity that is too low or too high may increase cross-talk between MRM channels. The Pulse Voltage is set to different values depending of the MS method being employed; the default values are shown in Table C-1. If the pulse voltage is set too low, cross-talk will increase and ions will slip over the pulse instead of “surfing” in front of it. Setting the pulse voltage to too high a value will result in a decrease in sensitivity.

Table C-1 Default Pulse Voltage Values

MS Mode	Pulse Voltage (V)	
	Collision Cell	Source
MS1	0 (T-Wave is "OFF")	2*
MS2	0 (T-Wave is "OFF")	2
Daughter	2	2
MRM	2	2
Parent	5	2
Neutral Loss/Gain	5	2
Survey	2	2

*** Note:** The source pulse voltage is only applied when using the optional MUX-technology Interface.