


Waters

MassLynx
Quantitation (New Version)

- Determines the concentration of specific analytes within a sample
- Can be done on data acquired through a variety of Acquisition Modes:
 - Multiple Reaction Monitoring (MRM)
 - Single Ion Recording (SIR)
 - Full Scan Acquisition
- QuanLynx and TargetLynx with an EPCAS System is Designed to Be a Part of a 21 CFR Part 11 Compliant Environment.

- In addition to unknown samples, a set of standards is also run to form a calibration curve.
- MassLynx analyzes the response of unknown samples and compares their response to that indicated by the calibration curve, then calculates the concentrations of the unknowns.



Steps in Creation of a Calibration Curve for Quantitation

- Integrate peaks in chromatograms
- In each chromatogram, determine the location of the peak relating to a specific compound
- Calculate response factor for the located peak
- Create a Calibration Curve for that compound

Project = Quantify.pro

Set of analyses on samples using a MS method that had:

- MRM of 3 channels
 - Internal Standard (294.1 > 64.0)
 - Analyte 1 - Parent Drug (288.1 > 58.0)
 - Analyte 2 - Metabolite (274.1 > 182.1)

- Used to account for Experimental Drift
- Can be Added at Various Points in the Analysis
 - - In the Original Sample
 - - Before Injection by the LC
- Response of Analyte in a Sample is:

$$\frac{\text{(Peak Area of Analyte)}}{\text{(Peak Area of I.S.) / (Conc of I.S.)}}$$

1. Enter Sample Types & Concentrations into Sample List
2. Determine Correct Integration Parameters for Chromatogram Peaks.
3. Create Quantification Method.
4. Process Samples.
5. Check Results – Adjust if Needed.
6. Print Out Results – Save Results on in Report File.

1. Set up Sample List

©2005 Waters Corporation

MassLynx - Quantify - Quantify.spl

File View Run Help

Shortcut Queue Status

Queue Is Empty

Spectrum Chromatogram Map Edit Samples

	File Name	File Text	MS File	Inlet File	Bottle	Inj	Sample Type	Conc A
1	ASSAY01	plasma blank	DEFAULT	DEFAULT	1	10	Blank	0
2	ASSAY02	0.2pg/ml std	DEFAULT	DEFAULT	2	10	Standard	0.2
3	ASSAY03	0.5pg/ml std	DEFAULT	DEFAULT	3	10	Standard	0.5
4	ASSAY04	0.75pg/ml std	DEFAULT	DEFAULT	4	10	Standard	0.75
5	ASSAY05	1pg/ml std						
6	ASSAY06	1pg/ml std						
7	ASSAY06	2pg/ml std						
8	ASSAY07	5pg/ml std						
9	ASSAY08	10pg/ml std						
10	ASSAY09	15pg/ml std						
11	ASSAY10	0.3pg/ml QC						
12	ASSAY11	2pg/ml QC						
13	ASSAY12	12pg/ml QC						
14	ASSAY13	Rat sample 01						
15	ASSAY14	Rat sample 02	DEFAULT	DEFAULT	14	10	Analyte	0
16	ASSAY15	Rat sample 03	DEFAULT	DEFAULT	15	10	Analyte	0
17	ASSAY16	Rat sample 04	DEFAULT	DEFAULT	16	10	Analyte	0
18	ASSAY17	Rat sample 05	DEFAULT	DEFAULT	17	10	Analyte	0
19	ASSAY18	Rat sample 06	DEFAULT	DEFAULT	18	10	Analyte	0
20	ASSAY19	Rat sample 07	DEFAULT	DEFAULT	19	10	Analyte	0

Ready Instrument Not Present 0:0 Only Error Shutdown Enabled

TargetLynx: Edit Method

BioLynx: Process Samples

QuantLynx: View Results

ChromaLynx: Edit Quan-Optimize Method

Run Quan-Optimize

NeoLynx: View Optimization Results

Sample List from Quantify.pro project.

QuanLynx selected from Shortcut Bar


1. Set up Sample List

©2005 Waters Corporation

- Standard Sample list plus two additional categories:
 - Sample Type
 - Concentration A (B, C, D.....)

	File Name	File Text	MS File	Inlet File	Bottle	Inject Volume	Sample Type	Conc A
1	ASSAY01	plasma blank	DEFAULT	DEFAULT	1	10.000	Blank	0
2	ASSAY02	0.2pg/ml std	DEFAULT	DEFAULT	2	10.000	Standard	0.2
3	ASSAY03	0.5pg/ml std	DEFAULT	DEFAULT	3	10.000	Standard	0.5
4	ASSAY04	0.75pg/ml std	DEFAULT	DEFAULT	4	10.000	Standard	0.75
5	ASSAY05	1pg/ml std	DEFAULT	DEFAULT	5	10.000	Standard	1
6	ASSAY06	2pg/ml std	DEFAULT	DEFAULT	6	10.000	Standard	2
7	ASSAY07	5pg/ml std	DEFAULT	DEFAULT	7	10.000	Standard	5
8	ASSAY08	10pg/ml std	DEFAULT	DEFAULT	8	10.000	Standard	10


1. Set up Sample List – Adding Extra Columns



File Text	Sample Type	Add	Net File	Bottle	Inj μ
plasma blank	Blank	Insert	ULT	1	10
0.2pg/ml std	Standard	Delete	ULT	2	10
0.5pg/ml std	Standard	Fill	ULT	3	10
0.75pg/ml std	Standard	Clear	ULT	4	10
1pg/ml std	Standard	Column	ULT	5	10
2pg/ml std	Standard	Format	ULT	6	10
5pg/ml std	Standard	Sort	ULT	10	10
10pg/ml std	Standard	Number of Samples...	ULT	10	10
15pg/ml std	Standard	Number of Injections...	ULT	10	10
0.3pg/ml QC	QC		ULT	10	10
2pg/ml QC	QC	2	DEFAULT	DEFAULT	10
12pg/ml QC	QC	12	DEFAULT	DEFAULT	10

Use the “**Samples / Format / Load**” menu item and load in a format that already has these fields.

Alternatively, you can ‘right click’ on the sample list and use the “**Customize Display**” item on the ‘pop-up’ menu and add these columns.




Blank -	Solvent or matrix, insures that system is clean and/or shows endogenous material in sample.
Standard -	Sample of a known concentration, used to form calibration curve.
Analyte -	Sample of unknown concentration.
QC -	Quality Control - Known concentrations, used to test the validity and accuracy of the calibration curve.

1. Specify Sample Types and Concentrations

©2005 Waters Corporation

- Pull Down menu within the sample list. Specify whether the sample is a Blank, Standard, Analyte or QC.



Alternatively, just type the first letter of the sample type (example for 'Analyte' type 'A') and hit enter.

Concentration A or (B, C...)

- The known concentrations of Standards or QC's must be entered into this column.

Conc A
0.2
0.5
0.75
1
2
5
10
15
0.3
2
12

2. Determine Correct Integration Parameters for Chromatogram Peaks

©2005 Waters Corporation

- Go to the Sample List and highlight a Standard in the middle of the concentration range.
- Click on the Chromatogram button.

The screenshot shows the MassLynx software interface. The title bar reads "MassLynx - Quantify - Quantify.spl". The menu bar includes "File", "View", "Run", and "Help". The toolbar contains icons for file operations and a "Shortcut" button. The main area is titled "Queue Is Empty" and features a navigation pane on the left with buttons for "Edit Method", "Process Samples", and "View Results". The "View Results" button is highlighted with a red arrow. The main table, titled "Samples", has columns for "File Name", "File Text", "MS File", "Inlet File", "Bottle", "Inj", "Sample Type", and "Conc A". The row for "ASSAY07" (5pg/ml std) is highlighted in black. The "Chromatogram" button in the top navigation bar is also highlighted with a red box and a red arrow.


	File Name	File Text	MS File	Inlet File	Bottle	Inj	Sample Type	Conc A
1	ASSAY01	plasma blank	DEFAULT	DEFAULT	1	10	Blank	0
2	ASSAY02	0.2pg/ml std	DEFAULT	DEFAULT	2	10	Standard	0.2
3	ASSAY03	0.5pg/ml std	DEFAULT	DEFAULT	3	10	Standard	0.5
4	ASSAY04	0.75pg/ml std	DEFAULT	DEFAULT	4	10	Standard	0.75
5	ASSAY05	1pg/ml std	DEFAULT	DEFAULT	5	10	Standard	1
6	ASSAY06	2pg/ml std	DEFAULT	DEFAULT	6	10	Standard	2
7	ASSAY07	5pg/ml std	DEFAULT	DEFAULT	7	10	Standard	5
8	ASSAY08	10pg/ml std	DEFAULT	DEFAULT	8	10	Standard	10
9	ASSAY09	15pg/ml std	DEFAULT	DEFAULT	9	10	Standard	15

Ready Instrument Not Present 0:0 Only Error Shutdown Enabled

2. Peak Integration-Display All Traces


- The TIC for the highlighted sample will be brought up.
- Click **Display, Mass** from the top of the chromatogram window.

click **Add Trace** and **Select All** to bring up all of the transitions.



2. Three Ion Chromatograms Should Now Be Shown

©2005 Waters



- Each individual transition will now be displayed.
- Delete the TIC to simplify the screen. (If the TIC is still displayed)

2. Setup Peak Integration-Noise

©2005 Waters Corporation


To setup the Integration use the (***Process, Integrate***) menu item to get the 'Integrate Chromatograms' dialog box. First determine the baseline noise by grabbing some noise (right click and drag) over a quieter area of the chromatogram .

(for example)

Peak to Peak Amplitude of Noise Will Be Filled In

Next Click on Smooth

Note Apex Track Peak Integration is Available




2. Setup Peak Integration-Smoothing

©2005 Waters Corporation

Continue setting up Integration Process: After Clicking on **Smooth**, Right click and drag over the peak at half height.


Window Size
Will Be
Filled In



Use Mean
Method When
Smoothing
Chromatograms

Remember the correct window size


2. Setup Peak Integration – Peak Detect




Without Apex
Peak Integration


Setups Baseline

Shoulder







If Apex Peak Integration is selected, 'noise' is handled in the peak detect setup. Noise and Peak width can be entered or you can have them calculated for you.



'Check' if you want Shoulders detected

Set starting and ending Baseline levels




Peak Threshold parameters can also be adjusted in the Method Editor. Specify Criteria to Discriminate Peaks from Noise.

For 'Relative' values, enter the percentage of the largest peak (base peak) a peak must exceed to be integrated.

For example, check 'Abs area' and enter 20% of the peak area for your lowest standard. Any peak with an area lower than this will be considered noise and not integrated.

- Click **OK**, the peak of interest will be integrated.
- Review the integration - is it acceptable? If not, repeat the integration with different parameters (noise, peak detect, thresholding) until satisfactory results are obtained.
- Once an acceptable integration is attained, you may want to test it on a low range standard and a high range standard to insure that parameters are adequate for the full range of response.




Here's an example of a well integrated peak and a poorly integrated peak.

2. Review Peak Integration- Example

©2005 Waters Corporation

Example of Peak Integration for Quantify Project

Note Smoothing Parameters for each transition



The screenshot displays the Waters MassLynx software interface. On the left, the 'Quantify' sidebar is visible, featuring a vertical menu with icons and labels for 'TargetLynx', 'BioLynx', 'QuanLynx', 'ChromaLynx', and 'NeoLynx'. The 'QuanLynx' section is highlighted, and the 'Edit Method' button is selected, indicated by a red arrow. The main window, titled 'Qmeth1.mdb - QuanLynx Method Editor', shows a 'Compound List' on the left and a 'Property' table on the right. The 'Compound List' contains three entries: '1: I. Std', '2: Parent', and '3: Metabolite'. The 'Property' table lists various method parameters and their values.

Property	Value
Compound Name	I. Std
Quantification Trace	294.10 > 64.00
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1.0000
Chromatogram mass window (PPM)	10.0000
Response Type	Internal (relative)
Internal Standard: 1	None
Internal Standard: 2	None
Internal Standard: 3	None
Internal Standard: 4	None
Internal Standard: 5	None
Internal Standard: 6	None

Example of QuanLynx Method Editor.

Click on 'Edit Method' to start method editor.

The different compounds in the assay are listed on the left.

Click on a compound in the list and the parameters that describe how to quantitate that compound are listed on the right.

Qmeth1_QL.mdb - QuanLynx Method Editor


File Edit View Compound Help

Compound List

- 1: I. Std
- 2: Parent
- 3: Metabolite

Property	Value
Compound Name	I. Std
Quantification Trace	294.10 > 64.00
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1.0000
Chromatogram mass window (PPM)	10.0000
Response Type	Internal (relative)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Fixed
Concentration of Standard	1.0000
<input type="checkbox"/> View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.2000

Use the buttons on the tool bar to decide which parameters for a compound you wish to view. 'Click' on the button on the left to display all of the parameters for a compound. 'Click' on one of the other buttons to display only a subset of the parameters.



All Params
(Display all
params)

Compound
Params
(e.g. Name,
m/z Trace)

Calculation
Factors
(Line fitting
params)

Integration
Params
(e.g.
smoothing,
peak detect)


Target Ion
Params
(confirma-
tory ions)

Calculation
Factors
(e.g. S/N)

For this example, use the 'File/New' menu item and create a new method.


'Click' on the 'Compound Properties' button to display the name and 'trace' fields.

To add a compound use the drop down menu or use the 'Add' button on the tool bar.



Compound Buttons

Import Add Delete




'Import' allows you to import info on a compound in from another method.


For this example, we are going to first enter quant parameters for the internal standard, next enter quant parameters for the parent drug and take care of the metabolite last.

5pg/ml std


ASSAY07 Sm (Mn, 2x2)



ASSAY07 Sm (Mn, 2x2)



ASSAY07 Sm (Mn, 2x2)



Hint. It will be easier if you:


©2005 Waters Corporation

Reduce the size of the Sample List window so this window along with the Method Editor dialog box occupies half of the screen and the chromatogram window occupies the other half of the screen.

Main MassLynx
Sample List Window

Method Editor
Dialog Box

Chromatogram
Window (pick a
standard from the
middle of the conc.
range and display
the channels from all
of the compounds.)



Size the windows so the Chromatogram window does overlap the other windows. (a little overlap is okay)

3. Build Quantitation Method Add First Compound to a New method

©2005 Waters Corporation

Click on the 'Compound/Add' menu item or the 'Add' button to add a compound to this method.

Example of new method after clicking on 'Add'.

The screenshot shows the 'Untitled - QuanLynx Method Editor' window. The 'Compound List' on the left contains one entry: '1: New Compound'. The 'Property' table on the right is as follows:

Property	Value
Compound Name	New Compound
Quantification Trace	
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	0.0200

Next type in the name for your compound (we will do the "Int Std" first in this example)

Example after entering name for our compound.

The screenshot shows the 'Untitled - QuanLynx Method Editor' window. The 'Compound List' on the left now contains one entry: '1: Int Std'. The 'Property' table on the right is the same as in the previous screenshot, but the 'Compound Name' value has been updated to 'Int Std'. A red arrow points from the text above to the 'Int Std' value in the table, which is also enclosed in a red box.

Property	Value
Compound Name	Int Std
Quantification Trace	
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	0.0200

Untitled - QuanLynx Method Editor

File Edit View Compound Help

Compound List
1: Int Std

Property	
Compound Name	
Quantification Trace	
Include Primary Trace in Response?	
Use absolute mass window?	
Chromatogram mass window (Da)	
Chromatogram mass window (PPM)	
Response Type	
Response Uses	
Acquisition Function Number	
Concentration of Standard: Level	
Concentration of Standard	
<input type="checkbox"/> View Retention Time Parameters	
Locate Peak Using	
Predicted Retention Time	
Predicted Relative Retention Time	
Relative Retention Time Reference	
Retention Time Window (mins) ±	
Lower Retention Time Tolerance (%)	
Upper Retention Time Tolerance (%)	
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest
<input type="checkbox"/> View Internal Standards	
Internal Standard: 1	None

Select Parameters to View

- Compound Name
- Quantification Trace
- Include Primary Trace in Response?
- Use absolute mass window?
- Chromatogram mass window (Da)
- Chromatogram mass window (PPM)
- Response Type
- Response Uses
- Acquisition Function Number
- Concentration of Standard: Level
- Concentration of Standard
- View Retention Time Parameters
- View Acceptance Flag Parameters

OK Cancel

To control the parameters that are displayed, 'Right Click' on the display. A window similar to the one at the right will appear and you can select which parameters you want displayed

©2005 Waters Corporation

'Right Click & Drag'

The screenshot shows the QuantLynx Method Editor interface. On the left, three chromatograms are displayed, each with a peak highlighted and its retention time and area noted. On the right, a 'Property' table is populated with values corresponding to the selected peak. Red arrows and boxes highlight the data transfer process.

Property	Value
Compound Name	Int Std
Quantification Trace	294.1 > 64
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	0.0200
Chromatogram mass window (PPM)	10.0000
Response Type	External (absolute)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Fixed
Concentration of Standard	0.0000
View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.7950
Predicted Relative Retention Time	0.0000
Relative Retention Time Reference	None
Retention Time Window (mins)	0.5500
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest
View Acceptance Flag Parameters	
View Internal Standards	
View Totals	

Chromatogram 1 (Top): 5pg/ml std, ASSAY07 Sm (Mn, 2x2), MRM of 3 Channels AP+, 2.79:844, 294.1 > 64, 4.51e3 Area, Int Std

Chromatogram 2 (Middle): ASSAY07 Sm (Mn, 2x2), MRM of 3 Channels AP+, 2.79:2415, 288.1 > 58, 1.21e4 Area

Chromatogram 3 (Bottom): ASSAY07 Sm (Mn, 2x2), MRM of 3 Channels AP+, 2.65:1302, 274.1 > 182.1

To enter parameters on the chromatogram trace for this compound, 'Right Click' and drag across the chromatographic peak for this compound. Trace, Acquisition Function, Ret Time and RT window info will be entered for you.

Compound Properties

- **Response Type**
–Set to “External (absolute)”
- **Concentration of Standard**
–Set Level to “Fixed”
– Set Concentration to “1”
- **Internal Standards**
–Set to “None” for the internal standard

Property	Value
Compound Name	Int Std
Quantification Trace	294.1 > 64
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1
Chromatogram mass window (PPM)	10.0000
Response Type	External (absolute)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Fixed
Concentration of Standard	1
<input type="checkbox"/> View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.7950
Predicted Relative Retention Time	0.0000
Relative Retention Time Reference	None
Retention Time Window (mins) ±	0.5500
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest
<input type="checkbox"/> View Acceptance Flag Parameters	
Maximum Blank Acceptance Response	<input checked="" type="checkbox"/> 0.0000
Maximum Concentration Limit	<input checked="" type="checkbox"/> 0.0000
Reporting Concentration Limit	<input checked="" type="checkbox"/> 0.0000
Minimum Recovery Level (%)	0.0000
Maximum Recovery Level (%)	100.0000
Flag Recovery Level:	<input checked="" type="checkbox"/> NO
<input type="checkbox"/> View Internal Standards	
Internal Standard: 1	None
Internal Standard: 2	None
Internal Standard: 3	None
Internal Standard: 4	None
Internal Standard: 5	None
Internal Standard: 6	None

Calibration Parameters

- Calibration Reference Compound**
 –Select the reference compound to match the current compound
- Polynomial Type**
 –Select “Average RF”
- Origin, Weighting Function, Axis Transformation**
 –These settings are not applied if the Polynomial Type is “Average RF” (which is the case for the internal standard)
- Concentration Units**
 –Type concentration unit of internal Standard (in this example it is pg/ml)
- Propagate Calibration Parameters**
 –When using an Internal Standard, the Propagate function is disabled
 –Click on the box until a red “X” appears. The Value should change to “No”

Property	Value
Compound Name	Int Std
Calibration Reference Compound	1: Int Std
Polynomial Type	Average RF
Origin	Exclude
Weighting Function	1/X
Axis Transformation	None
Concentration Units	pg/ml
User RF Value	<input checked="" type="checkbox"/> 0.0000
Propagate Calibration Parameters	<input checked="" type="checkbox"/> NO

Integration Parameters

- Smoothing**
 - To apply smoothing, click on the box until a green “Check” appears. The Value should change to “Yes”
 - Select the appropriate Smoothing Method, and enter the Smoothing Iterations and Width determined from the Peak Integration performed in the Chromatogram Window
- Threshold Parameters**
 - To apply threshold, click on the box next to the appropriate threshold parameter until a green “Check” appears. Enter the threshold value to be applied
- Propagate Integration Parameters**
 - Under most conditions, the integration parameters can be propagated for all compounds
 - To Propagate, click on the box until a green “Check” appears. The Value should change to “Yes”

Property	Value
Compound Name	Int Std
Smoothing Enabled?	<input checked="" type="checkbox"/> YES
<input type="checkbox"/> View Smooth Parameters	
Smoothing Method	Mean
Smoothing Iterations	1
Smoothing Width	2
Apex Track Enabled?	<input checked="" type="checkbox"/> NO
<input type="checkbox"/> View Peak Detect Parameters	
<input type="checkbox"/> Apex Track Parameters	
Peak-to-Peak Baseline Noise	<input checked="" type="checkbox"/> 10.0000
Peak Width at 5% Height	<input checked="" type="checkbox"/> 30.00
Baseline Start Threshold %	0.05
Baseline End Threshold %	0.05
Detect Shoulder Peaks?	<input checked="" type="checkbox"/> NO
<input type="checkbox"/> Standard Peak Detection Parameters	
Peak-to-peak noise amplitude	0.0000
Automatic Noise Measurement	<input checked="" type="checkbox"/> YES
Balance	30
Splitting	90
Detect Shoulder Peaks?	<input checked="" type="checkbox"/> NO
Detect Shoulder Peaks Threshold	30
Reduce Tail	100
Reduce Height	5
<input type="checkbox"/> View Threshold Parameters	
Threshold Relative Height	<input checked="" type="checkbox"/> 1.5000
Threshold Absolute Height	<input checked="" type="checkbox"/> 0.0000
Threshold Relative Area	<input checked="" type="checkbox"/> 2.0000
Threshold Absolute Area	<input checked="" type="checkbox"/> 0.0000
Integration Window Extent	0.0000
Propagate Integration Parameters?	<input checked="" type="checkbox"/> YES

Integration Parameters

- Peak Detect Parameters**

– To enable Apex Track, click on the box until a green “Check” appears. The Value should change to “Yes” (In this example Apex Track is disabled)

– Apex Track Parameters: Typical setting are shown. These parameters are applied only when Apex Track is enabled

– Standard Peak Detection Parameters: Typical settings are shown. These parameters are applied when Apex Track is disabled

Property	Value
Compound Name	Int Std
Smoothing Enabled?	<input checked="" type="checkbox"/> YES
<input type="checkbox"/> View Smooth Parameters	
Smoothing Method	Mean
Smoothing Iterations	1
Smoothing Width	2
Apex Track Enabled?	<input type="checkbox"/> NO
<input type="checkbox"/> View Peak Detect Parameters	
<input type="checkbox"/> Apex Track Parameters	
Peak-to-Peak Baseline Noise	<input checked="" type="checkbox"/> 10.0000
Peak Width at 5% Height	<input checked="" type="checkbox"/> 30.00
Baseline Start Threshold %	0.05
Baseline End Threshold %	0.05
Detect Shoulder Peaks?	<input type="checkbox"/> NO
<input type="checkbox"/> Standard Peak Detection Parameters	
Peak-to-peak noise amplitude	0.0000
Automatic Noise Measurement	<input checked="" type="checkbox"/> YES
Balance	30
Splitting	90
Detect Shoulder Peaks?	<input type="checkbox"/> NO
Detect Shoulder Peaks Threshold	30
Reduce Tail	100
Reduce Height	5
<input type="checkbox"/> View Threshold Parameters	
Threshold Relative Height	<input type="checkbox"/> 1.5000
Threshold Absolute Height	<input type="checkbox"/> 0.0000
Threshold Relative Area	<input type="checkbox"/> 2.0000
Threshold Absolute Area	<input type="checkbox"/> 0.0000
Integration Window Extent	0.0000
Propagate Integration Parameters?	<input checked="" type="checkbox"/> YES

- **Secondary Ion Parameters**

–Secondary Ion Trace: If a secondary ion is being used, then the transition is entered (in this example there is no Secondary Ion Trace)

–Use trace in response calculation: If enabled, then the area of the secondary trace is added to the area of the primary trace when calculating the response. To enable, click on the box until the green “Check” appears. The value should change to “Yes”

–Secondary Ion Ratio: Enter the ratio of the secondary ion to the primary ion. This can be calculated from the chromatogram window

Target Ion Parameters

Property	Value
Compound Name	I. Std
<input checked="" type="checkbox"/> View Secondary Ion Parameters	
Secondary Ion Trace	
Use trace in response calculation?	<input checked="" type="checkbox"/> NO
Secondary Ion Ratio	0.0000
Secondary Ion Ratio Tolerance (%)	0.0000
Secondary Ion Must Exist?	<input checked="" type="checkbox"/> NO
Secondary Ion Must Pass Ratio?	<input checked="" type="checkbox"/> NO

3. Quant Method Editor – Enter Calculation Factors

Calculation Factors

- **Detection Limit Parameters**

- Signal-to-noise method: If RMS is used, then Noise calculation factor should be “3”; If Peak-to-peak is used, then Noise calculation factor should be “1”

- Noise window: Specify portion in chromatogram from which noise will be calculated. If both fields are set to “0.0000”, then the noise window will be determined automatically by the software

- Detection/Quantification Limit Factor: User-defined parameters. Typical Values are shown

Property	Value
Compound Name	Int Std
+ View Toxic Equivalence Factors	
+ View Mole Ratio	
- Detection Limit Parameters	
Signal-to-noise method	RMS
Noise calculation factor	3.0000
Noise window start (min)	0.0000
Noise window end (min)	0.0000
Measure peak signal level from	
Detection Limit Factor	3
Quantitation Limit Factor	8
Propagate Detection Limit Settings?	<input checked="" type="checkbox"/> YES
Use EMPC?	<input type="checkbox"/> NO
User Peak Factor	0.0000

This can be:

- A single decimal number (m/z) for mass chromatograms (from SIR or Full Scan (continuum or centroid))
- Two decimal numbers separated by a ">" for an MRM function e.g. 609.2 > 195.1
- 'TIC' for total ion current chromatograms
- 'BPI' for base peak intensity chromatograms
- An1, An2, An3, or An4 for analog data
- The wavelength for DAD data.
- Ch1, Ch2 etc for SIR data to use one quantify method with multiple SIR functions. Where Ch1 is the first mass in the list, Ch2 is the second etc.

View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.7950
Predicted Relative Retention Time	0.0000
Relative Retention Time Reference	None
Retention Time Window (mins) ±	0.5500
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest

Peak Location and Time Window
Parameters can be entered from the keyboard if needed.

- ***Peak Location: Retention Time (RT) and Time Window***

–Parameters were entered during the ‘right click and drag’ over the peak. ‘RT’ is center of a time interval that the peak must appear in to be associated with this compound. ‘Time Window’ is the width of this interval. So for this example, the peak must appear at 2.7950 ± 0.55 min.

For 'Concentration of Standards' select 'Conc A' (or 'Conc B', 'Conc C', etc.) from the drop down menu for the column in the sample list that has the concentration values for this compound.

If this compound is an internal standard you may select 'Fixed' (usually done for Int Std's). If you enter 'Fixed', you can enter the concentration of the Int Std or since it is the same concentration in all of the samples, simply enter 1.000 (usually done for Int Std's).

Property	Value
Compound Name	Int Std
Quantification Trace	294.1 > 64
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1
Chromatogram mass window (PPM)	10.0000
Response Type	External (absolute)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Fixed
Concentration of Standard	1

[-] View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.7950
Predicted Relative Retention Time	0.0000
Relative Retention Time Reference	None
Retention Time Window (mins) ±	0.5500
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest 

- **Peak Selection**

- If more than one peak is detected in the 'Time Window', this designates which peak to choose.
- The peak **Nearest** to the entered RT,
- the **Largest** peak,
- the **First** peak,
- the **Last** peak
- or **Totals** (sum up all of the peaks).

If your compound always has the largest peak in the chromatogram, select '**Largest**'.


- This entire process now needs to be repeated for the two other compounds (for the Quantify.pro example).
- Things that may differ between compounds:
 - Transition (Quantify Trace)
 - Name
 - Integration Parameters
 - Internal Reference (Select Internal Standard if Used)
 - Concentration of Standards
 - Retention Time
 - Time Window
 - Response Type in General Parameters Window
 - Polynomial Type in the General Parameters Window

- Click **Add New Compound Icon** to enter the compound into the Compound List (Shown is New Compound entered into the list)
- Note: Information from the previous compound has been propagated to the new compound. Check all parameters to ensure proper quantification of the new compound


Compound List	Property	Value
1: Int Std	Compound Name	New Compound
2: New Compound	Quantification Trace	294.1 > 64
	Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
	Use absolute mass window?	<input checked="" type="checkbox"/> YES
	Chromatogram mass window (Da)	1.0000
	Chromatogram mass window (PPM)	10.0000
	Response Type	External (absolute)
	Response Uses	Area
	Acquisition Function Number	One
	Concentration of Standard: Level	Fixed
	Concentration of Standard	1.0000
	<input type="checkbox"/> View Retention Time Parameters	
	Locate Peak Using	Retention Time
	Predicted Retention Time	2.7950
	Predicted Relative Retention Time	0.0000
	Relative Retention Time Reference	None
	Retention Time Window (mins) ±	0.2300
	Lower Retention Time Tolerance (%)	0.0000
	Upper Retention Time Tolerance (%)	0.0000
	Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
	Locate Peak Selection	Nearest
	<input checked="" type="checkbox"/> View Acceptance Flag Parameters	
	<input checked="" type="checkbox"/> View Internal Standards	
	<input checked="" type="checkbox"/> View Totals	

5pg/ml std


ASSAY07 Sm (Mn, 2x2) MRM of 3 Channels AP+
 294.1 > 64
 4.51e3
 Area



ASSAY07 Sm (Mn, 2x2) MRM of 3 Channels AP+
 288.1 > 58
 1.21e4
 Area



ASSAY07 Sm (Mn, 2x2) MRM of 3 Channels AP+
 274.1 > 182.1
 6.72e3
 Area



Property	Value
Compound Name	Parent
Quantification Trace	288.1 > 58
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1.0000
Chromatogram mass window (PPM)	10.0000
Response Type	Internal (relative)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Conc A
Concentration of Standard	1.0000
<input type="checkbox"/> View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.7950
Predicted Relative Retention Time	0.0000
Relative Retention Time Reference	None
Retention Time Window (mins) ±	0.5500
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest
<input type="checkbox"/> View Acceptance Flag Parameters	
<input type="checkbox"/> View Internal Standards	
<input type="checkbox"/> View Totals	

- Right Click and Drag on the chromatogram for the new compound trace
- Rename the new compound

3. Quant Method Editor – Enter Compound Properties (Analyte)

Compound Properties

- **Response Type**
–Set to “Internal (relative)”
- **Concentration of Standard**
–Set Level to “Conc A (or B,C,...)”
- **Internal Standards**
–Select the compound from the ‘drop down’ list that is the internal standard (In this case, Int Std).

Property	Value
Compound Name	Parent
Quantification Trace	288.1 > 58
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1.0000
Chromatogram mass window (PPM)	10.0000
Response Type	Internal (relative)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Conc A
Concentration of Standard	1.0000
<input type="checkbox"/> View Retention Time Parameters	
Locate Peak Using	Retention Time
Predicted Retention Time	2.7950
Predicted Relative Retention Time	0.0000
Relative Retention Time Reference	None
Retention Time Window (mins) ±	0.5500
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest
<input type="checkbox"/> View Acceptance Flag Parameters	
Maximum Blank Acceptance Response	<input checked="" type="checkbox"/> 0.0000
Maximum Concentration Limit	<input checked="" type="checkbox"/> 0.0000
Reporting Concentration Limit	<input checked="" type="checkbox"/> 0.0000
Minimum Recovery Level (%)	0.0000
Maximum Recovery Level (%)	100.0000
Flag Recovery Level?	<input checked="" type="checkbox"/> NO
<input type="checkbox"/> View Internal Standards	
Internal Standard: 1	1: Int Std
Internal Standard: 2	None
Internal Standard: 3	None
Internal Standard: 4	None
Internal Standard: 5	None
Internal Standard: 6	None

3. Quant Method Editor – Enter Calibration Parameters (Analyte)

Calibration Parameters

- **Polynomial Type**
 - For Analytes, “Linear” is the typical polynomial type. Other options include “Quadratic, Cubic, Quartic”
- **Concentration Units**
 - Type concentration unit of internal Standard (in this example it is pg/ml)
- **Propagate Calibration Parameters**
 - When using an Internal Standard, the Propagate function is disabled
 - Click on the box until a red “X” appears. The Value should change to “No”

Property	Value
Compound Name	Parent
Calibration Reference Compound	2: Parent
Polynomial Type	Linear
Origin	Exclude
Weighting Function	1/X
Axis Transformation	None
Concentration Units	pg/ml
User RF Value	<input checked="" type="checkbox"/> 0.0000
Propagate Calibration Parameters	<input checked="" type="checkbox"/> NO


3. Quant Method Editor – Enter Calibration Parameters (Analyte)

Calibration Parameters

- Origin, Weighting Function, Axis Transformation**
 - Origin: Typically the origin is excluded as a point of the calibration curve
 - Weighting Function: A weighting factor is appropriate for a calibration range that is greater than an order of magnitude. For this example, a weighting of “1/X” will be used
 - Axis Transformation: When the weighting function is used, the Axis Transformation function is not applied.

Property	Value
Compound Name	Parent
Calibration Reference Compound	2: Parent
Polynomial Type	Linear
Origin	Exclude
Weighting Function	1/X
Axis Transformation	None
Concentration Units	pg/ml
User RF Value	<input checked="" type="checkbox"/> 0.0000
Propagate Calibration Parameters	<input checked="" type="checkbox"/> NO

©2005 Waters Corporation



Property	Value
Compound Name	Metabolite
Quantification Trace	274.1 > 182.1
Include Primary Trace in Response?	<input checked="" type="checkbox"/> YES
Use absolute mass window?	<input checked="" type="checkbox"/> YES
Chromatogram mass window (Da)	1.0000
Chromatogram mass window (PPM)	10.0000
Response Type	Internal (relative)
Response Uses	Area
Acquisition Function Number	One
Concentration of Standard: Level	Conc A
Concentration of Standard	1.0000
<input type="checkbox"/> View Retention Time Parameters	
Locate Peak Using	Relative Retention Time
Predicted Retention Time	2.65
Predicted Relative Retention Time	.950
Relative Retention Time Reference	2: Parent
Retention Time Window (mins) ±	0.3200
Lower Retention Time Tolerance (%)	0.0000
Upper Retention Time Tolerance (%)	0.0000
Flag RT Tolerance?	<input checked="" type="checkbox"/> NO
Locate Peak Selection	Nearest
<input type="checkbox"/> View Acceptance Flag Parameters	
<input type="checkbox"/> View Internal Standards	
<input type="checkbox"/> View Totals	

R.T. of Metabolite Relative to Parent
 $(2.65) / (2.79) = 0.950$

Example of Use of
 Relative Retention Time

©2005 Waters Corporation


Minor Points About the Method Editor.

User RF Value

If no calibration curve is available, then divide the peak area (response) by this factor to calculate the conc.'s.


User Peak Factor

All concentrations calculated will be multiplied by this factor.



The screenshot shows the 'User RF Value' property in the Quant Method Editor. A red box highlights the 'User RF Value' property, and a red arrow points from the text 'User RF Value' to this property. The 'User RF Value' is set to 0.0000 with a red 'X' icon next to it. The 'Propagate Calibration Parameters' property is also highlighted with a red box and has a red 'X' icon next to it, set to 'NO'.

Property	Value
Compound Name	Int Std
Calibration Reference Compound	1: Int Std
Polynomial Type	Average RF
Origin	Exclude
Weighting Function	1/X
Axis Transformation	None
Concentration Units	pg/ml
User RF Value	<input checked="" type="checkbox"/> 0.0000
Propagate Calibration Parameters	<input checked="" type="checkbox"/> NO



The screenshot shows the 'User Peak Factor' property in the Quant Method Editor. A red box highlights the 'User Peak Factor' property, and a red arrow points from the text 'User Peak Factor' to this property. The 'User Peak Factor' is set to 0.0000. The 'Propagate Detection Limit Settings?' property is checked and set to 'YES'.

Property	Value
Compound Name	Int Std
+ View Toxic Equivalence Factors	
+ View Mole Ratio	
- Detection Limit Parameters	
Signal-to-noise method	RMS
Noise calculation factor	3.0000
Noise window start (min)	0.0000
Noise window end (min)	0.0000
Measure peak signal level from	
Detection Limit Factor	3
Quantitation Limit Factor	8
Propagate Detection Limit Settings?	<input checked="" type="checkbox"/> YES
Use EMPC?	<input checked="" type="checkbox"/> NO
User Peak Factor	0.0000

- Once the entire method is built, it's time to process the samples.
- Highlight the samples to quantitate. If the entire sample list is to be processed, click on the upper left box to activate the entire sample list
- Click on *QuanLynx, Process Samples*.

MassLynx - Quantify - Quantify.spl

File View Run Help

Queue Is Empty

QuanLynx ?

Edit Method

Process Samples

View Results

	File Name	File Text	Sample Type	Conc A
1	ASSAY01	plasma blank	Blank	0
2	ASSAY02	0.2pg/ml std	Standard	0.2
3	ASSAY03	0.5pg/ml std	Standard	0.5
4	ASSAY04	0.75pg/ml std	Standard	0.75
5	ASSAY05	1pg/ml std	Standard	1
6	ASSAY06	2pg/ml std	Standard	2
7	ASSAY07	5pg/ml std	Standard	5
8	ASSAY08	10pg/ml std	Standard	10
9	ASSAY09	15pg/ml std	Standard	15
10	ASSAY10	0.3pg/ml QC	QC	0.3
11	ASSAY11	2pg/ml QC	QC	2

Highlight the samples you wish to quantitate and then select 'Process Samples'

4. Processing Samples with QuanLynx

©2005 Waters Corporation

This window appears to confirm the specifics prior to processing samples. Double-check that the:

- 1) Designated samples are correct
- 2) That it is using the correct method.

For complete quantitation set this window to:

- 1) Integrate the chromatograms
- 2) Create a calibration curve
- 3) Calc the concentrations in each sample

The screenshot shows the 'Create QuanLynx Dataset' dialog box. On the left, there is a list of options with checkboxes and icons:

- Update Method Times
- Integrate Samples
- Calibrate Standards
- Quantify Samples
- Print Quantify Reports
- Export Results to LIMS

On the right, there are input fields and buttons:

- Project:** C:\MassLynxBackup_5\Quantify.PRD
- Quantify:**
 - From Sample: 1
 - To Sample: 39
 - Method: Example_Method (with a Browse button)
 - Curve: My_Great_Curve (with a Browse button)
- LIMS Export:**
 - File: Lims (with a Browse button)

At the bottom right are 'OK' and 'Cancel' buttons. Red arrows point from the text above to the 'From Sample' and 'To Sample' fields, and from the 'Integrate Samples' checkbox to the 'Integrate Samples' text in the list.

After samples are done processing, a “**Quantify**” box will appear on the lower tool bar. Double click to bring up the results.

OR

Using the main toolbar, click on the **View Results** button.

©2005 Waters Corporation

From the View drop down menu you set this window to show (or not show):


Information Bar (name of compound)

Summary Table

Chromatograms

Calibration Curves

'Double Click' on a Sample in the list to see its chromatograms




5. Some Features of the Menu of the QuanLynx Quantification Results Viewer


©2005 Waters Corporation

Step Through
the Different
Samples

Step Through
the Different
Compounds




Turn On or Off
the Slide Show



Experimental Record Can Be
Displayed Using Menu Item

5. Manually Adjusting Integrations

©2005 Waters Corporation




Click on the Baseline & 'Grab' the End Point using the Pointer and Left Mouse Button & Move the End Point to the desired spot.

A faint line will show the position of the original baseline. Reports will now show that this baseline was manually adjusted.

5. Manually Adjusting Integrations

©2005 Waters Corporation




'Right Click' in the window and the 'pop up menu' shown to left will appear.

Select 'Save Peak Mod' from this menu to 'accept' and 'save' the changed baseline.

Go to another chromatogram to keep the original baseline.

5. Editing Calibration Curves



To remove a point from the calibration, 'Right Click' on the point. The dialog box shown below should appear.

Show Chromatograms

Exclude

You can 'Left Click' and 'rubber band' zoom in on a region of the curve to get a closer look at small regions of the curve. Use the  button at the top to 'Unexpand'.


Note you can also 'Left Click' and 'rubber band' zoom in on a regions of the chromatograms.

'Left Click' on 'Exclude' and the point will be removed from the curve. Reverse the process to put the point back into the curve.

©2005 Waters Corporation

If you have the residuals displayed in your calibration window, you can also remove a point from the calibration by 'Right Clicking' on the 'bad' point.

'Right Click' on the point you wish to exclude here.



5. Right Click on the Display and Select 'Display Options' to Customize the Display:

©2005 Waters Corporation

For Example the Chromatogram Display Can Be Changed:


Use 'Summary' Page to Control Slide Show

Chromatogram Display Adjustments

Toggle Showing Internal Standard Chromatogram

Y Axis


Peak Annotation




X Axis

5. 'Right Click' On the Summary Table and ...


1) Click on 'Change Column Order' and which columns are displayed in the table can be changed.




2) Click on 'Edit Column Properties' and the properties of which ever column you 'right clicked' on can be changed.



From the 'File' Menu, Select the 'Report Format' Item:



For Example the Compound & Sample Summaries Can be Formatted Differently From How They Appear on the Screen




Customized Display and Report Formats Can Be Saved For Later Use

- Printing Reports (**File, Print Report**). Besides a full report, results from a set range of samples can be printed.
- Screen and Report Format. A customized format can be saved in a *.fmt file for later use.
- The Quantify Method used with a report can be changed using (**Edit, Quantify Method**)
- Editing of Calibration Curve (**Edit, Calibration Curve**) allows excluding of specific data points. ('Right Click' on a point in a Calibration Curve and select 'Exclude Point').
- Reprocessing samples after editing Quantify Method (**Process, Calculate**)

- Everything is in One File
- This file can be viewed and reports printed at a later date without reprocessing data
- This File Will Contain:
 - Compound and Sample Summaries
 - Calibration Curves
 - Chromatograms
 - Experimental Record for Each Analysis Run
 - Quantitation Method


Linear Least Squares Line Fitting and 1/X Weighting Factors



Conc	Meas.	Fitted	Diff	Diff ²
2	200	200	0	0
3	300	300	0	0
4	400	400	0	0
5	540	500	40	1600

Sum of
Diff² = 1600


Conc= Std Concentration
 Meas=Measured Peak Area
 Fitted=Peak Area from Fitted
 Line at given Conc
 Diff= Meas – Fitted
 Diff2= Diff Squared



Conc	Meas.	Fitted	Diff	Diff ²
2	200	192	8	64
3	300	304	- 4	16
4	400	416	- 16	256
5	540	528	12	144

Sum of
Diff² = 480

Conc= Std Concentration
 Meas=Measured Peak Area
 Fitted=Peak Area from Fitted
 Line at given Conc
 Diff= Meas – Fitted
 Diff2= Diff Squared



Using
'Possible'
Fitted Line

Conc	Meas.	Fitted	Diff	Diff ²
2	200	200	0	0
3	300	300	0	0
4	400	400	0	0
5	540	500	40	1600


Sum of
Diff² = 1600


Using
Linear
Least
Squares
Fitted Line

Conc	Meas.	Fitted	Diff	Diff ²
2	200	192	8	64
3	300	304	- 4	16
4	400	416	- 16	256
5	540	528	12	144


Sum of
Diff² = 480

Conc	Peak Area
1	9
3	30
10	100
30	300
100	1050









Diff = 1 Diff ² = 1






Diff = 50 Diff ² = 2500

Closer Look at Linear Least Squares Line Fitted with No Weighting









Nom Conc	Peak Area	Line Fitted w/ No Weighting		Line Fitted w/ 1/X Weighting	
		Calc Conc	Rel Error	Calc Conc	Rel Error
1	10	1.45	45.0%	1.12	12.4%
3	30	3.35	11.7%	3.05	1.6%
10	100	10.00	0.0%	9.78	-2.2%
30	300	29.01	-3.3%	29.01	-3.3%
100	1000	95.53	-4.5%	96.33	-3.7%