GE Power & Water
Water & Process Technologies
Analytical Instruments

Americas: T 800 255 6964, 303 444 2009 | geai@ge.com

Europe, Middle East, Africa: T 44 (0) 161 864 6800 | geai@europe@ge.com

Asia Pacific: geai.asia@ge.com

GE Power & Water
Water & Process Technologies
Analytical Instruments

QUICK START GUIDE

GE Autosampler

for use with Sievers M-Series Laboratory and Portable Analyzers

©2016 General Electric Company. All rights reserved.

DQS 77700-01_MUL Rev B

QUICK START GUIDE

GE Autosampler

for use with Sievers* M-Series Portable and Laboratory TOC Analyzers

Contents:

Specification Summary p. 3

Installation

English	p. 4-5
Português	p. 6-7
中文 (Chinese)	p. 8-9
Español	p. 10-11
日本語 (Japanese)	p. 12-13
Italiano	p. 14-15
Deutsch.....	p. 16-17
Français	p. 18-19

Warnings

English	p. 20
Português	p. 20
中文 (Chinese)	p. 21
Español	p. 21
日本語 (Japanese)	p. 22
Italiano	p. 22
Deutsch	p. 23
Français	p. 23

IMPORTANT!

For more detailed instructions, you can download the **Sievers M-Series Operation and Maintenance Manual** from our website. Go to www.geinstruments.com/library/manuals, and then make your selection.

* Trademark of General Electric Company; may be registered in one or more countries.

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Specification Summary

GE Autosampler for Use with M-Series Laboratory and Portable Analyzers	
Operating Specifications	
Dimensions (Base Model)	52.17 cm height x 29.2 cm width x 53.3 cm depth (20.5 in. height x 11.5 in. width x 21.0 in. depth)
Dimensions (with Rinse Station)	52.1 cm height x 31.1 cm width x 53.3 cm depth (20.5 in. height x 12.25 in. width x 21.0 in depth.)
Net Weight	Approximately 15 kg (33 lb)
Power Supply	100-240 VAC, 160 VA, 50-60 Hz
Fuse Rating	Replace with same type and size fuse, T1.6A
Inputs/Outputs	USB — Used for connecting to an M-Series Analyzer. RS-232 — Used for connecting to a previous generation Analyzer.
Temperature Range (Non-Operating)	5-45°C, 20-95% RH (41-113°F, 20-95% RH)
Temperature Range (Operating)	10 °C - 40 °C, 8 - 95% RH (50 °F - 104 °F, 8 - 95% RH)
Racks – Maximum #	Sample racks -3 (plus 1 Emergency rack)
Certification	CE

¹ The power supply is a recognized cRU component for Canada and the U.S.

QUICK START GUIDE

ENGLISH

Installation Instructions

STEP 1: Unpack and Inspect the GE Autosampler

STEP 2: Complete the Identification Records

Enter the serial number of the GE Autosampler on page 3 of the *Sievers M-Series Operation and Maintenance Manual*.

STEP 3: Select a Location


Place the GE Autosampler next to an installed Analyzer on a clean, unobstructed surface that can accommodate the dimensions and weight detailed in the Specifications Summary. For proper heat dissipation, ensure that 16 cm (6.3 in) is available at the rear and on both sides of the Analyzer.

Keep away from direct sunlight. For proper operation and to avoid measurement errors, ensure that the environment temperature is no greater than 40° C (104° F) or lower than 10° C (50° F).

STEP 4: Set up the GE Autosampler

Install the rinse station tubing (optional), vial racks, and needle assembly.

Install the Rinse Station Tubing (Optional)

Routing for Rinse Station Tubing

Intall the Racks


Insert each vial rack into the base of the GE Autosampler, fitting the feet (located on the bottom of the vial rack) into the corresponding slots.

Installed Vial Racks

Note: All racks in a single protocol must be either 17 mL, 35 mL, 40 mL, or 60 mL.

Install the Needle Assembly

Needle Arm and Needle Assembly

1. Using a flat-blade screwdriver (**turn counter-clockwise ONLY**), slightly lower the needle arm to make room for using the wrench in Step 4.
2. Unscrew the PEEK nut at the top of the needle assembly and remove the protective cap from needle tip.

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WARNINGS ENGLISH

For continued protection against fire hazard, replace fuse with same type and rating.

- This is a Safety Class 1 product. It must be attached to a grounded power source. If this instrument is used in a manner not specified by GE Analytical Instruments USA, the protection provided by the instrument may be impaired.
- Always stop TOC measurements/protocols before turning off or unplugging the GE Autosampler.
- Observe the following guidelines when using the GE Autosampler in conjunction with flammable liquids. To reduce the risk of fire, the GE Autosampler shall not be used to sample flammable liquids with a flash point of less than ambient temperature + 25 °C. Other than the needle, the GE Autosampler should have no parts in contact with the fluid. Ensure that no parts, in or near the fluid surface are capable of igniting a flammable liquid.
- Keep clear of the probe arm during operation. Although it is designed to stop upon detecting an obstruction, the probe arm can move without warning during normal operation and at start-up.

This symbol on the instrument indicates that the product does not contain any restricted substances included in China RoHS I.

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3. Position the needle into the arm, lowering the tip through the lower mount block, and pushing the top of the assembly upward and into the upper mount block.
4. Secure with the nut (placed over the washer), and tighten with 1/2" open-end wrench. Be careful NOT to over-tighten the nut.
5. Screw the PEEK nut (with sample tubing attached) back into the threaded mount of the needle assembly until it is finger tight. DO NOT over-tighten the nut.
6. Using a flat-blade screwdriver, raise the arm back into place.


STEP 5: Route the Sample Tubing

1. Insert the anchor guide into the hole on the top side of the GE Autosampler.

Anchor Guide


2. Route the tubing through the small eyelet on the Autosampler arm.


Small Eyelet on Needle Arm

3. Connect the Valco fitting to the end of the tubing.
4. Route the tubing through the anchor guide.

STEP 6: Connect the Communication/Power Cables

Connecting Directly to a Network

Connecting to a Computer

1. Connect the GE Autosampler to the Analyzer using the provided USB-to-USB cable.
2. Connect the Analyzer to a network using a straight-thru Ethernet cable or directly to a computer using a USB cable. Insert the power cord into the GE Autosampler and then connect to a grounded power outlet.

STEP 7: Install DataPro2 Software

Refer to the *GE Autosampler Installation Guide* for instructions. (Download from www.geinstruments.com/library/manuals).

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