GE Power & Water
Water & Process Technologies
Analytical Instruments

QUICK START GUIDE

M-Series TOC Analyzers

M9, M9^e, and M5310 C (On-Line, Laboratory, and Portable)

QUICK START GUIDE

M-Series TOC Analyzers

M9, M9^e, and M5310 C (On-Line, Laboratory, and Portable)

Contents:

Specification Summary p. 3

Installation

English p. 4-5
Deutsch p. 6-7
Français..... p. 8-9
Italiano..... p. 10-11
中文 (Chinese) p. 12-13
日本語 (Japanese) p. 14-15
Español p. 16-17
Português p. 18-19

Warnings

English p. 20
Deutsch..... p. 21
Français..... p. 22
Italiano..... p. 23
中文 (Chinese) p. 24
日本語 (Japanese)..... p. 25
Español p. 26
Português p. 27

IMPORTANT!

For more detailed instructions, you can download the **Operation and Maintenance Manual** for your Analyzer model from our website. Go to www.geinstruments.com/library/manuals, and make your selection.

* Trademark of General Electric Company; may be registered in one or more countries.

Specification Summary

	M9 and M9 ^e Laboratory Analyzer	M9 and M9 ^e On-Line Analyzer	M9 and M9 ^e Portable Analyzer
Operating Specifications			
Range	0.03 ppb to 50 ppm		
Precision	<1% RSD		
Accuracy	± 2% or ± 0.5 ppb, whichever is greater		
Sample Type	Autosampler or discrete grab sample	On-line continuous or discrete grab sample	On-line continuous, Autosampler, or discrete grab sample
Display Readout	3 significant digits		
Calibration	Typically stable for 12 months		
Analysis Time	2 minutes (4 seconds with the optional Turbo mode)		
Sample Temperature ¹	5-95 °C (41-203 °F)		
Ambient Temperature	5-40 °C (41-104 °F)		
Sample Pressure	n/a	100 psig	
On-Line Flow Rate	n/a	>50 mL/min (for on-line mode)	
Instrument Sample Flow Rate	0.5 mL/min		
Analyzer Specifications			
Outputs	USB device port (1), USB host ports (3); Modbus TCP/IP	4-20 mA outputs (3); alarm outputs (4); binary input (1); USB device port (1), USB host ports (2);Modbus TCP/IP	
Display	7" WVGA 800x480 pixel, Color LCD w/touch-screen		
Power	100 – 240 V~, 50 – 60 Hz, 100 VA		
Fuses	Replace with same type and size fuse: T 1.6 A 250 VAC Fuse (Slow Blow), size 5 x 20 mm appliance inlet		
Dimensions	H: 42.2 cm (16.6 in); W: 24.6 cm (9.7 in); D: 40.0 cm (15.8 in)	H: 54.9 cm (21.6 in); W: 45.0 cm (17.7 in); D: 26.5 cm (10.4 in)	H: 39.5 cm (15.4); W: 22.9 cm (9.0 in); D: 46.4 cm (18.3 in)
Weight	9.4 kg (20.6 lb)	15.8 kg (34.9 lb)	9.4 kg (20.8 lb)
Enclosure Rating	n/a	IP-45	IP-21
Safety Certifications	ETL CE		
Sample Conductivity ² (Optional)			
Raw Conductivity Range	0.01 to 2,000 µS/cm	n/a	0.01 to 2,000 µS/cm
Conductivity Accuracy	±0.005 µS/cm or ±1%, whichever is greater	n/a	±0.005 µS/cm or ±1%, whichever is greater
Conductivity Precision	<0.25% RSD	n/a	<0.25% RSD
Environment			
Maximum Relative Humidity	0 – 95%, non-condensing		
Maximum Altitude	3,000 m (9,800 ft)		
Pollution Degree	2		

	M5310 C Laboratory Analyzer	M5310 C On-Line Analyzer	M5310 C Portable Analyzer
Operating Specifications			
Range	4 ppb to 50 ppm		
Precision	<1% RSD		
Accuracy	± 2% or ± 0.5 ppb, whichever is greater		
Sample Type	Autosampler or discrete grab sample and TOC removal grab	On-line continuous or discrete grab sample, timed on-line, TOC removal on-line, TOC removal grab	On-line continuous, Autosampler, or discrete grab sample, timed on-line, TOC removal on-line, TOC removal grab
Display Readout	3 significant digits		
Calibration	Typically stable for 12 months		
Analysis Time	2 minutes		
Sample Temperature	5-60 °C (41-140 °F)		
Ambient Temperature	5-40 °C (41-104 °F)		
Sample Pressure	n/a	100 psig	
On-Line Flow Rate	n/a	>50 mL/min (for on-line mode)	
Instrument Sample Flow Rate	0.5 mL/min		
Analyzer Specifications			
On-Line Inlet(s)	n/a	One stream, or two-stream inlet (option)	One stream
Outputs	USB device port (1), USB host ports (3); Modbus TCP/IP	4-20 mA outputs (3); alarm outputs (4); binary input (1); USB device port (1), USB host ports (2); Modbus TCP/IP	
Display	7" WVGA 800x480 pixel, Color LCD w/ touch-screen		
Power	100 – 240 V~, 50 – 60 Hz, 100 VA		
Fuses	Replace with same type and size fuse: T 1.6 A 250 VAC Fuse (Slow Blow), size 5 x 20 mm appliance inlet		
Dimensions	H: 42.2 cm (16.6 in.); W: 24.6 cm (9.7 in); D: 40.0 cm (15.8 in)	H: 54.9 cm (21.6 in); W: 45.0 cm (17.7 in); D: 26.5 cm (10.4 in)	H: 39.5 cm (15.4); W: 22.9 cm (9.0 in); D: 46.4 cm (18.3 in)
Weight	9.7 kg (21.4 lb)	16.2 kg (35.6 lb)	9.8 kg (21.6 lb)
Enclosure Rating	n/a	IP-45	IP-21
Safety Certifications	ETL CE		
Environment			
Maximum Relative Humidity	0 – 95%, non-condensing		
Maximum Altitude	3,000 m (9,800 ft)		
Pollution Degree	2		

¹ If the sample temperature is above 60°C (140°F), the optional PVDF iOS is required.

² Sample conductivity is available in Autosampler or discrete grab sample mode only with GE Dual-Use Conductivity TOC vials.

QUICK START GUIDE

ENGLISH

Installation Instructions

Note: For more detailed instructions, refer to the Analyzer's Operation and Maintenance Manual, downloadable from www.geinstruments.com (Library > Manuals).

Step 1: Unpack and Inspect the Analyzer

Step 2: Complete the Identification Records

Complete page 3 of the *Operation and Maintenance Manual*.

Step 3: Select a Location

Place the Analyzer away from direct sunlight and avoid temperatures greater than 40 °C (104 °F) or lower than 10 °C (50 °F).

Mount the On-Line TOC Analyzer on a wall or a support stand with the screen at eye level (or set it on a workbench that can support its weight). Select a stand or bench, and mounting bolts, to support four times the weight of the instrument [total = 64.8 kg (142.40 lb)]. Allow a minimum of 5 cm (2 in.) clearance between the back of the Analyzer and the wall and 30.5 cm (12.0 in.) on the sides, top, and bottom.

Place the *Portable* or *Laboratory* TOC Analyzer on a clean surface that can accommodate the Analyzer's dimensions and weight (See the [Specification Summary](#)).

Step 4: Install the Reagent Cartridges


Reagent Cartridge

Note: Installation instructions are for the type of reagent cartridge shown here. If installing the other type with the red-and-green button valve, refer to the Analyzer's Operation & Maintenance Manual.


1. Ensure power is OFF and exercise ESD protection.
2. Open the Analyzer door (*On-Line*) or remove the right side panel (*Laboratory* and *Portable*).
Stabilize the Analyzer with one hand on the top panel (*Laboratory Analyzer*) or handle (*Portable Analyzer*). Grasp the edge of the right side panel from the back of the Analyzer. Pull the entire panel away from the Analyzer. (The panel is attached to the Analyzer by magnets.)
3. Remove the reagent tray (*Portable* and *Laboratory* models) or remove the reagent tray cover (*On-Line* model).
4. Pull the yellow vent sticker off the bottom of the *Acid* cartridge.
5. With the cartridge upright, remove the containment tape from the inlet nozzle and remove the plug.
6. Locate the *Acid* supply line, attach it to the top of the *Acid* reagent cartridge, and tighten finger tight.
7. Place the *Acid* cartridge in the back compartment of the reagent tray, ensuring that the nozzle-end points out through the cutout.
8. Repeat steps 4-7 for the *Oxidizer* cartridge.

9. For *Portable* or *Lab* models, place the reagent tray onto the unit floor. For *On-Line* models, carefully replace the cover.

Laboratory* Reagent Tray


*Portable reagent tray position is similar to Laboratory shown.


On-Line Reagent Tray

Step 5: Install the DI Water System

1. Assemble the water bottle and insert the fill fitting into the **Prime** port.
2. Fill the DI reservoir with deionized water. (Do not overfill. Ensure that the reservoir float is just below the top of the window).


Filling the DI Loop

Reservoir Float

3. Install the resin cartridge into the **PRIME** position to ready for priming in Step 11.

PRIME Position

Step 6: Install I/O Control Wiring (*On-Line* and *Portable* Only)

Refer a qualified electrician to the Analyzer's *Operation and Maintenance Manual* for installing I/O wiring, alarms, and remote start binary input.

Step 7: Install Power

Connect the power cord to *Laboratory* and *Portable* Analyzers. Or for wiring AC conduit to the On-Line Analyzer, see the Analyzer's *Operation and Maintenance Manual*.

Do NOT turn the power on until instructed in Step 11.

Step 8: Install a Printer (Optional)

Install a printer (Bixolon/Samsung SRP-275, Citizen CD-5500A, or Epson TM-7220D).


1. Connect the printer cable's USB host connector into one of the **USB host** ports on the Analyzer.
2. Connect the cable's USB device connector to the printer.
3. Configure the printer to the Analyzer later in Step 12.

Step 9: Install the iOS System

(When connecting the Portable or On-Line Analyzer for online measurements.)

Note: If using a GE Autosampler, DO NOT proceed with this section. Instead, refer to the GE Autosampler Installation Guide.

1. Ensure that the power to the Analyzer is off and disconnected.
2. Verify that the water source is disabled.
3. Connect the 1/4" PFA inlet tubing (w/ in-line filter) into the sample inlet **(A)** on the iOS System. Secure nut with 9/16" open-end wrench. Do NOT over-tighten.
4. For a two-stream option (M5310 C On-Line Only), connect the other inlet tubing into the second sample inlet.
5. Slide one end of the waste tubing over the waste outlet fitting **(B)** and position the tubing for a gravity drain.
6. Ensure that the stainless steel tubing is connected to the **Sample Inlet** and establish a water flow to the iOS System.
7. Measure the flow rate and adjust to between 50-300 mL/min using the flow rate adjustment screw **(C)**.

iOS

Step 10: Install the Waste Outlet

(When configuring the Laboratory for grab or Autosampler use.)

1. Using a 7/16" open-end wrench, disconnect the nut that connects to the **Waste** port.
2. Attach the 1/8" OD plastic waste tubing to the **Waste** port, then use the 7/16" open-end wrench to tighten securely.
3. Position the hose for a gravity drain of the waste water.


Step 11: Prime the DI Water Pump

Note: The DI water reservoir was filled in Step 5 and the DI pump is now ready for priming.

1. Ensure that the resin cartridge is in the **Prime** position (Step 5).
2. Press the power switch to *On* and let the Analyzer run 180 min.
3. After 180 minutes, on the *Home* screen, shut down the system and turn off power to the Analyzer.

Warning: To avoid the possibility of corrupting the Analyzer's operating system, always **FIRST** shut down the system on the Analyzer's *Home* screen before pushing the Analyzer's power switch to the *Off* position.

1. View the reservoir water level and verify that the float is positioned **JUST BELOW** the top of the reservoir. (The color of the resin in the cartridge should be black.)
2. If the water level is now lower due to the circulation process, fill the DI water reservoir as needed. There is no need to prime again.
3. Install the DI resin cartridge into the **RUN** position. Ensure that the ends are secure.
4. Close and latch the Analyzer's door (*On-Line* model) or replace the panel (*Laboratory* and *Portable* models).

RUN Position

Step 12: Configure Basic Analyzer Settings

Before placing the Analyzer into service, record the installation of the reagents on the *Consumables* tab of the *Maintenance* screen. In the appropriate *Acid* or *Oxidizer* section, press *RENEW*.

Use the Analyzer screen to configure various basic settings, including optional security (*Password Protection* or *DataGuard*) using the touchscreen display. Most of these settings will not need to change unless moving the Analyzer or reconfiguring the operational environment.

Step 13: Flush the Analyzer

Use these steps to flush the reagent syringes.

To also rinse the Analyzer's sample fluidics (optional for low TOC applications), refer to the *Analyzer's Operation and Maintenance Manual* for instructions.

To flush the reagent syringes

1. Connect to either a continuous DI water source or a large flask (1000-2000 mL) of DI water.
2. Power on the Analyzer.
3. On the *Maintenance* screen, select the *Advanced* tab.
4. Press **Syringe Flush** and follow the on-screen prompts.
5. Shut down the system and turn the power switch to *Off*.

Step 14: Install a GE Autosampler

(Laboratory and Portable TOC Analyzers Only)

If planning to use a GE Autosampler with the Analyzer, install it now. Refer to the *GE Autosampler Installation Guide* for complete installation instructions. (Download from www.geinstruments.com (Library>Manuals).

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Faça o download em www.geinstruments.com (Library > Manuals).

WARNINGS

This symbol on the instrument indicates that the user should refer to the manual for operating instructions.

- Shipping boxes containing an Analyzer may weigh over 40 lb —necessitating a two-person lift technique when handling. Such boxes are labeled with this warning, supporting OSHA standards intended to protect employees from injury.
- If this instrument is used in a manner not specified by GE Analytical Instruments, the safety protection provided by the instrument may be impaired.
- To avoid injury and to ensure accurate analysis results, ensure that the Analyzer door (On-Line TOC Analyzer) is closed and latched or side panels (Laboratory and Portable TOC Analyzers) are on and in place during normal operation.
- The Analyzer requires an A/C mains supply of 100-240 volts and requires the correct power cord for safe operation. Use only the supplied power cord with the Analyzer, which meets this A/C mains supply requirement.
- For safety reasons and to avoid the possibility of corrupting the database, ALWAYS follow this procedure when shutting down the Analyzer. First stop analysis AND stop the operating system (using the Shut down button in the Analyzer's firmware) BEFORE pressing the Analyzer's power switch to Off or unplugging the power cord.
- Any operation requiring access to the inside of the Analyzer, including installation of maintenance items, could result in injury. To avoid potentially dangerous electrical shock, take the following precautions. BEFORE opening the Analyzer, stop any analysis, shut down the operating system (using the Shut down button in the Analyzer's firmware) and turn off power and disconnect from the power supply.
- Electrical conduit and output and alarm connections should be installed by a qualified electrician. ESD protection should also be exercised.
- This is a Class A product. In a domestic environment, this product may cause electromagnetic interference. Take adequate measures to correct the interference.
- Hazardous reagents (ammonium persulfate and phosphoric acid) are used in the Analyzer. Before installing a reagent, read the corresponding Material Safety Data Sheet (MSDS) for proper handling precautions and spill or leak procedures. The MSDS are contained in the pouch located on the outside of the reagent shipping box.
- Installation of reagents and most Analyzer maintenance requires access to the inside of the Analyzer. To avoid potentially dangerous shock do this. BEFORE opening the Analyzer, stop any analysis and shut down the operating system (using the Shut Down button in the Analyzer's firmware), then turn off power and disconnect from the power supply.
- To avoid exposure to the chemical reagents, wear acid-resistant gloves and safety goggles.
- Reagent cartridges are for single-use only. Do NOT refill. Any refilling or reusing of reagent cartridges will void all Analyzer and parts warranties and nullify any performance claims.
- After removing the plug on the BLOW MOLD type reagent cartridge, keep the cartridge in an upright position before attaching and detaching the reagent supply line to avoid chemical from dripping out and avoiding exposure to the chemical reagent. Save the plug after removing to later reinsert into the nozzle inlet after removing the cartridge from the Analyzer during a cartridge replacement process.
- To protect against accidental exposure to ultra-violet radiation, do not operate the UV lamp outside of its protective housing.
- (On-Line and Portable TOC Analyzers) Water in the iOS System may be hot. Before inserting a vial into the iOS System after operating in On-Line mode, slide the door open and wait 30 seconds to allow sample to completely drain. Inserting a vial before draining can result in hot water spray projecting upward out of the iOS System.
- (Laboratory and Portable TOC Analyzers) For continued protection against fire hazard, replace fused with the same type and rating.

(On-Line TOC Analyzer) This symbol indicates the protective earth terminal (ground) for the Analyzer.

The iOS System and vial ports contain sharp needles designed to pierce the septa of sample vials. Do NOT put fingers or inappropriate materials into the iOS System or vial port.

When servicing parts inside the Analyzer, ensure that power to the Analyzer is OFF, and keep hands clear of the reagent syringe assemblies. The syringes are controlled by moving parts that can pinch skin.

HAZARDOUS WASTE DISPOSAL

Hazardous reagents (ammonium persulfate and phosphoric acid) are used in the Analyzer. The waste stream from the instrument is acidic and must be disposed of properly. Consult your federal, state, and local government regulations.

- The UV lamp and the display screen contain mercury and may be considered hazardous material in your local area. Dispose of these items in accordance with federal, state, or local government regulations.
- In the case of any broken or damaged UV lamp, handle the remains in accordance with your organization's toxic waste handling procedure and dispose of them in accordance with federal, state, or local government regulations.

This symbol indicates the European Union Directive 2002/96/EC for waste electrical and electronic equipment (WEEE) requires disposing the Analyzer separately from standard waste.

OPERATIONAL CAUTIONS

- Ensure that the DI water Reservoir is full, particularly when running samples with high TOC or high salt concentrations. After running high TOC or salt samples, always rinse the Analyzer by running low-TOC DI water through a sample cycle.
- To avoid false TOC readings and possible damage to the Analyzer, always ensure that the sample inlet is open and the DI water reservoir is filled before starting analysis.
- Operation of the Analyzer without the in-line filter on the sample inlet line will damage the Analyzer and void the warranty. To avoid damaging the Analyzer, install the filter and replace the filter element, as needed.
- To avoid damage to the DI water pump, ensure that the DI water reservoir is filled BEFORE turning On the Analyzer.
- To avoid false TOC readings and possible damage to the Analyzer, always ensure that the sample is flowing through the iOS System and the DI water reservoir is filled before starting analysis.
- For performance within specifications on ozonated water systems, an Ozone Destruct Kit must be purchased from GE Analytical Instruments and installed according to instructions.
- When cleaning the Analyzer, always apply liquid to the cloth directly.

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©2014 General Electric Company. All rights reserved.